

APPLICATION NOTE

**ECLIPSE Ji Digital Inverted Microscope
AX Confocal Microscope System
General Analysis 3 Image analysis Software module**

3D image analysis reveals drug-induced cellular senescence and increased nucleolus volume: Efficient imaging methods to maximize data acquisition



High content analysis incorporates multiple measurements such as cell morphology and fluorescence intensity of specifically-labelled cellular features to assess biological activity in response to set of conditions, and is a powerful approach for applications like drug compound screening. Wide-field microscopy offers relatively high imaging throughput, which is ideal for screening drug candidate compounds, while confocal laser scanning microscopy enables three-dimensional visualization of structures and volumetric quantification to obtain deeper insight into compound responses. Using both wide-field and confocal laser scanning microscopy, drug discovery and life science research can be advanced more efficiently. In this application note, we introduce an imaging method that efficiently utilizes both 2D and 3D analysis of wide-field and confocal images to reveal how drug-induced cellular senescence impacts the nucleolus.

Keywords: cellular senescence, high content analysis, wide-field microscopy, confocal microscopy, 3D analysis, volume measurement, nucleolus, DNA damage, γ H2AX, granule analysis, drug discovery

Sample preparation and experimental results

2D analysis (Fig. 1): A549 cells were treated with Doxorubicin (DOX) at concentrations of 50, 100, 200, and 400 nM to induce senescence. The senescent A549 cells were seeded in a 96-well plate, and nuclei, nucleoli, and γ H2AX granules were fluorescently stained. First, we identified compound concentrations that showed changes by 2D analysis of nucleoli and γ H2AX using a Wide-field microscope.

3D analysis (Fig. 2): Next, we used a confocal microscope to acquire images of DOX concentrations (0, 100, 400 nM) and ran 3D analysis of nuclei and nucleoli. 2D analysis showed that the area of the nucleoli were similar at different DOX concentrations (Fig.1 (E)). However, 3D analysis revealed that the volume of the nucleoli increased at higher compound concentrations (Fig. 2(E)). This clearly indicates that 3D analysis of nucleoli volume is a critical analytical parameter in the study of cellular senescence.




Fast imaging with a large field of view and high resolution to capture tiny granules of γ H2AX

Wide-field

ECLIPSE Ji built-in camera and 20x lambda D objective combine to achieve a resolution per pixel of **0.316 $\mu\text{m}/\text{pixel}$**

Even small granular structures of less than 1 μm in size can be resolved with a low-magnification 20x objective, allowing the capture of a large number of cells in a single shot with the potential for quantitative statistical analysis.



Objective: PLAN APO λ D 20x (NA0.8)
Built-in camera, Wide-field (D-LEDI2)

2D & 3D high-content analysis to maximize data acquisition efficiency

Wide-field

Confocal



Wide-field: Study populations

- ✓ Drug Screening
- ✓ Drug concentration considerations



Confocal: Explore details

- ✓ Gain subcellular insight
- ✓ Make new discoveries
- ✓ Conduct 3D analysis

Imaging Modality	Nucleus	Nucleolus	γ H2AX	Features
Wide-field (2D)	Area	Fluorescence Intensity	<ul style="list-style-type: none"> • Number • Area • Fluorescence Intensity 	<ul style="list-style-type: none"> ✓ Population overview ✓ Fast image acquisition
Confocal (3D)	Volume	Volume		<ul style="list-style-type: none"> ✓ 3D image ✓ Volume info

Table 1: Analytical parameters and optimal imaging method features for studying cellular senescence

STEP 1 Wide-field: Screening

- ✓ Rapidly capture images
- ✓ Low-magnification objective captures many cells for quantitative statistical analysis
- ✓ Good for drug screening
- ✓ Lighter data volume and easy analysis

2D: 20x Air objective

Low Magnification Statistical analysis

Fast

Constructing EDF images from Z-stack images for 2D analysis

[DOX]: 0, 50, 100, 200, 400 nM
Experimental groups: 5 conditions
2 images/well, 15 wells
Total: 30 images

STEP 2 Confocal: Details

- ✓ Three-dimensional visualization
- ✓ 3D analysis
- ✓ Volume information

3D: 40x Water immersion objective

High magnification 3D analysis

3D analysis by constructing 3D images from Z-stack images

[DOX]: 0, 100, 400 nM
Experimental groups: 3 conditions
3 images/well, 3 wells
Total: 9 images

Detection region	Fluorescence label	Ex/Em (nm)
Nucleus	DAPI	345/455
Nucleolus	Nucleolus Bright Red	537/605
γ H2AX (DNA damage)	DNA Damage Detection Kit - γ H2AX - Deep Red	646/668
Light source and detector		
Epi-fluorescence (D-LEDI2), Ji Built-in camera (FOV 25mm)		
Objective (N.A.)		
CFI Plan Apochromat Lambda D 20X (N.A. 0.8)		
Magnification	Field of view (FOV), Resolution	
20X	0.88 x 0.88 mm / image 2800 x 2800 pixel	
# of points	Z-stacks	
2 points/well	1.825 μm x 3 steps (Range: 3.75 μm)	
# of wells	Drug (Doxorubicin) concentrations	
15 wells	0, 50, 100, 200, 400 nM	
# of images	Image processing	
30 images	Constructing EDF images from Z-stack images	

Table 2: Detection area, fluorescent labels, and conditions for image acquisition using a wide-field microscope for Fig 1.

Detection region	Fluorescence label	Ex/Em (nm)
Nucleus	DAPI	345/455
Nucleolus	Nucleolus Bright Red	537/605
Light source and detector		
Confocal Laser (AX), DUX-VB 4Ch Detector Unit (FOV 25mm)		
Objective (N.A.)		
CFI Apo LWD Lambda S 40XC WI (N.A. 1.15)		
Magnification	Field of view (FOV), Resolution, Others	
40X	0.44 x 0.44 mm / image 1024 x 1024 pixel Galvano, Averaging: 2	
# of points	Z-stacks	
3 points/well	0.32 μm x 45-70 steps (Range: 14 - 22 μm)	
# of wells	Drug (Doxorubicin) concentrations	
3 wells	0, 100, 400 nM	
# of images	Image processing	
9 images	Constructing a 3D image from Z-stack images	

Table 3: Detection area, fluorescent labels, and conditions for image acquisition using a confocal laser-scanning microscope for Fig 2.

Fast, accurate imaging for new discoveries

In this experiment, we used a wide-field microscope to capture a large number of images in a short time, allowing us to study the effects of a drug on variably-treated populations across the entire plate. A single, fully in-focus EDF image was constructed from three Z-stack images, and the nuclei, nucleoli, and γ H2AX granules were binarized and measured. γ H2AX granules appeared with maximum count, area, and fluorescence intensity at 200 nM DOX treatment. On the other hand, the number of nucleoli and their fluorescence intensities were maximal with 100 and 400 nM Dox treatment, respectively. The area of nucleoli showed no difference with drug concentration above 50 nM. To then study subcellular details, confocal images of A549 cells at three different drug concentrations were acquired. The 3D confocal images were binarized and the volumes of the nuclei and nucleoli were measured. The results of the 3D analysis revealed that high concentrations of Doxorubicin increased nucleolus volume.

Summary

- ✓ Fast, 2D statistical analysis with wide-field microscopy
- ✓ Three-dimensional structural examination and 3D analysis with confocal laser scanning microscopy
- ✓ Volumetric information leads to new discoveries
- ✓ Efficient imaging was combined 2D and 3D analysis
- ✓ Nucleolus volume is a useful analytical parameter in the study of cellular senescence
- ✓ Combined imaging modalities allow accurate detection of drug effects
- ✓ Seamless 2D & 3D analysis from drug discovery to pharmacological effects with a single ECLIPSE Ji

Acknowledgments

Sample provided by Dojindo Laboratories, Co. Ltd. We would like to express our sincere gratitude to everyone at Dojindo Laboratories for their cooperation in establishing the experimental conditions for senescence induction and the staining conditions protocol optimized for imaging.

Sample Preparation Protocol

Please refer to the Nikon ECLIPSE Ji Application Note for the A549 cell senescence induction and sample preparation protocol "[Label-free quantitative analysis of cellular senescence and high-content imaging of granules](https://www.microscope.healthcare.nikon.com/ja_JP/resources/application-notes/label-free-quantitative-analysis-of-cellular-senescence-and-high-content-imaging-of-granules)"

https://www.microscope.healthcare.nikon.com/ja_JP/resources/application-notes/label-free-quantitative-analysis-of-cellular-senescence-and-high-content-imaging-of-granules

Cell staining protocol

■ DNA Damage Detection Kit - γ H2AX - Deep Red

* This is slightly different from the method described in the G267 instruction manual.

- 1) Remove the cell culture supernatant, add 4% PFA solution (100 μ l/well) to each well, and incubate at room temperature for 3 minutes.
- 2) Remove the PFA solution and wash the cells three times with PBS (100 μ l/well).
- 3) add 0.1% Triton X-100/PBS solution (100 μ l/well) and incubate at room temperature for 30 minutes.
- 4) Remove the supernatant and wash the cells twice with PBS (100 μ l/well).
- 5) Add Blocking Solution (100 μ l/well) and incubate for 20 minutes at room temperature.
- 6) Remove the supernatant and wash the cells twice with PBS (100 μ l/well).
- 7) Add γ H2AX staining solution (100 μ l/well) and incubate overnight at 4 degrees.
- 8) Remove the supernatant and wash the cells twice with PBS (100 μ l/well).
- 9) Add Secondary antibody staining solution (100 μ l/well) and incubate for 1 hour at room temperature.
- 10) Remove the supernatant and wash the cells twice with PBS (100 μ l/well).

■ Nucleolus Bright Red

- 1) After γ H2AX staining, remove the supernatant and add a mixture of adjusted Nucleolus Bright working solution and DAPI (100 μ l/well) and incubate at room temperature for 5 minutes. * Please refer to the instruction manual of N512. Nucleolus Bright Red:1000x dilution, DAPI:1000x dilution https://www.dojindo.co.jp/manual/N511_N512/
- 2) The supernatant is removed and the cells are washed twice with PBS (100 μ l/well).

Reagents, Product name	Product number	Supplier
-Cellstain®- DAPI solution	D523	Dojindo Laboratories
DNA Damage Detection Kit - γ H2AX - Deep Red	G267	Dojindo Laboratories
Nucleolus Bright Red	N512	Dojindo Laboratories

References

Abigail Buchwalter, *et al.*, Nucleolar expansion and elevated protein translation in premature aging. *Nature Communications* volume 8, Article number: 328 (2017)

Alba Corman, *et al.*, Targeting the nucleolus as a therapeutic strategy in human disease. *Trends Biochem Sci.* Volume 48, Issue 3, March 2023, Pages 274-287


Product information

ECLIPSE Ji Digital Inverted Microscope with AX Confocal Microscope System

The ECLIPSE Ji is capable of fast, high-resolution image acquisition using its built-in camera and epi-fluorescence light source. This digital inverted microscope is optimized for 2D high content analysis, with the ability to upgrade to confocal when required. From drug screening using fast wide-field imaging to understanding three-dimensional structures and 3D analysis using confocal images, you can obtain seamless results with this single unit.



3D mask image Mesh display



General Analysis 3 image analysis software module

By simply combining analysis blocks, nucleus and nucleolus regions can be easily binarized and measured, allowing flexible image analysis.

The "ConnectCells^{3D}" function allows you to easily perform 3D analysis in steps similar to a 2D analysis.

