
Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. August 2023 ©2023 NIKON CORPORATION

Monitor images are simulated.

TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.

Business Unit website



NIKON CORPORATION

Company names and product names appearing in this brochure are their registered trademarks or trademarks.

Shinagawa Intercity Tower C, 2-15-3, Konan, Minato-ku, Tokyo 108-6290, Japan 471, Nagaodai-cho, Sakae-ku, Yokohama, Kanagawa 244-8533 Japan phone: +81-3-6433-3705 fax: +81-3-6433-3785

https://www.healthcare.nikon.com/en/

ISO 14001 Certified for NIKON CORPORATION

Nikon Instruments Inc.

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A. phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only) fax: +1-631-547-0299A

https://www.microscope.healthcare.nikon.com/

Nikon Europe B.V.

Stroombaan 14, 1181 VX Amstelveen, The Netherlands phone: +31-20-7099-000

https://www.microscope.healthcare.nikon.com/en_EU/

Nikon Precision (Shanghai) Co., Ltd.

CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060 (Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026 (Guangzhou branch) phone: +86-2-3882-0551 fax: +86-2-3882-0580 https://www.microscope.healthcare.nikon.com/zh_CN/

Nikon Canada Inc.

CANADA phone: +1-905-625-9910 fax: +1-905-602-9953

Nikon France, Succursale de Nikon Europe B.V. FRANCE phone: +33-1-4516-4516

Nikon Deutschland, Zweigniederlassung der

Nikon Europe B.V.

GERMANY phone: +49-211-9414-888 Nikon Italy, Branch of Nikon Europe B.V.

ITALY phone: +39-055-300-9601

Nikon Europe B.V., Amstelveen, Zweigniederlassung

Schweiz (Egg/ZH)

SWITZERLAND phone: +41-43-277-2867

Nikon UK, Branch of Nikon Europe B.V. UNITED KINGDOM phone: +44-208-247-1717

Nikon Österreich, Zweigniederlassung der Nikon Europe B.V.

AUSTRIA phone: +43-1-972-6111

Nikon Singapore Pte. Ltd.

SINGAPORE phone: +65-6559-3651 fax: +65-6559-3668

Nikon Australia Pty Ltd AUSTRALIA phone: +61-2-8767-6900

Nikon Instruments Korea Co., Ltd.

KOREA phone: +82-2-6288-1900 fax: +82-2-555-4415

Nikon India Private Limited

INDIA phone: +91-124-4688-500

Printed in Japan (2308)D Code No. 2CE-MQER-2 This brochure is printed on recycled paper made from 40% used material.



DIGITAL IMAGING MICROSCOPE ECLIPSE Ui

AREAL SOLUTION FOR DIGITAL PATHOLOGY

Live on-screen observation.

Nikon's new ECLIPSE Ui Digital Upright Microscope provides accurate microscopy-based pathology imaging. View and share high quality images in real time with easy to use software for a simpler workflow.

VIEW

Display high quality images while remaining simple to operate

FAST

Immediate response for quicker workflow.

USABILITY

User-friendly to promote operational efficiency.

Sample images are monitored in real time. The quality is backed by Nikon's renowned imaging technology - clear color reproducibility without negative influence from ambient lights. Operator eye fatigue is greatly reduced as the need to look through conventional eyepieces is eliminated.

The system is operational in 2.5 seconds after loading a sample. Digital sample images can be observed live, plus magnification changes and XY movements can be quickly adjusted. It is also equipped with macro-imaging function and other sample-oriented applications.

The GUI (graphical user interface) is intended for easy identification and for efficient observation tasks. The functions needed to observe the sample images are arranged in an operator-friendly and efficient manner.

DAILY SUPPORT

Functions to support multiple use cases and applications.

The system is provided with three modes: routine specimen observation tasks, research* and education*, and data sharing. Users can select their preferred imaging quality and speed. Automatic bar code linkage from slide to image ensures sample control.

* Not for use in diagnostic procedures.



DIGITAL IMAGING MICROSCOPE

ECLIPSE Ui




VIEW

Display high quality images while remaining simple to operate.

New category medical equipment for turning pathological samples digital.

The internal PC provides all the necessary functions and applications.



2 Images on the monitor screen for easy observation.


It is no longer necessary to sit for hours looking through microscope eyepieces. The images are shared onscreen, suitable for two or more people to discuss the samples.



Time-proven optical-mechanical performance of Nikon microscopes.

The CFI Plan Fluor series objectives are backed by Nikon's superior optical technology. Their transmittance and numerical aperture (resolving power) are acclaimed for their high standard.





Immediate response for quicker workflow.

Operator-friendly

A prepared slide with sample can be placed with one hand. Onscreen control for magnification changes, X-Y movement and focusing makes for efficient performance.









Samples are loaded at the touch of a button and displayed on the monitor in 2.5 seconds.

The macro-imaging function is available for overviewing samples.

A sample is captured for a macro-image together with the target location. They are displayed along with the stage. Using the overview, the region in question can be observed. That site can be marked on the spot and the preset monitoring location can be recalled with just one click. The macro-image displayed may be saved as required.



Equipped with various support functions.

The system is designed to meet various use cases, such as thick sample observation and recording, and successive observation.

Z-stack function*

The Z-stack function enables to observing and recording images that are thick, undulating or disperse in the Z axis.

*Not for use in diagnostic procedures.

Stepwise transfer

For convenient successive observation, the stage can be shifted at a constant rate. This step-by-step function has been developed with less afterimage in mind. Six speed settings can be preset for suitable scrolling.

STEP STEP STEP STEP

Starting the system

- Press the Power button.
- Log in.
- Select an operating mode.

Preparations

- Place a prepared slide with sample on the holder.
- Press the Sample Load/Unload button.

- · Check the sample ID information.
- Look at the macro-image to check the full view.

Checking the entire image

- Send the macro-image to the pathological system.
- Specify the area to observe on the macro-image.

Observing live

- Focus adjustment.
- Exposure adjustment.
- Refer to pp. 7-8 for details.

Observation complete

- Press the Sample Load/Unload button.
- Take out the prepared slide.

Shutting down the system


• Press the Power button.

05


USABILITY To promote operational efficiency.


The system is equipped with the GUI for easy identification and for efficient observation tasks.


The functions needed to observe sample images are arranged in a user-friendly way. Simultaneously captured micro- and macro-images are displayed live.











Magnifications changeable with the touch of a button

Digital zoom-in/out available.

Easy Focus function

Equipped with the auto-focus function. Focus adjustment also possible with the onscreen button AF or the mouse wheel.

Brightness adjustment with the slide bar

Exposure and aperture adjustable by dragging the slide bars onscreen to the right or left.

Color adjustment

Shades and contrasts changeable as required.

Annotations added*

Areas of interest can be marked as well as point-to-point measurements made in the displayed image.

*Not for use in diagnostic procedures.



DAILY SUPPORT Functions for use cases, applications of images and sample control.

Mode selection to meet use cases.

Routine

In this mode, cytoscreening is supported. Sample images displayed live on the monitor screen are used for routine pathological image observation and successive observation. Images can also be transferred to pathological systems



Research*

Data (sample images, observation spots, etc.) to share or discuss are saved in an external storage*2. This data can be utilized for relevant studies and education.

- *1: Not for use in diagnostic procedures.
- *2: Separately sold.

ECLIPSE Ui Display

Remote

This mode allows remote user (users on contract) in remote locations to operate the system in real time. These users can also observe images*3.

*3: Separate contract must be concluded for using the Remote mode. For communication environment, contact us.



Compact size and a well though-out design

The space-saving body measures 422 mm in height, 233 mm in width and 427 mm in depth. Ambient light does not affect images. The microscope is readily set up.






5 Different resolutions and frame rates to meet applications.

There are two types of image output: LIVE for immediate observation and eveluation, and High-Quality for saving and storage.

LIVE (Observation mode) 1080 x 1080, above 30 fps

High-Quality (Capture mode) 3712 x 3712, above 1 fps



3712 x 3712, above 1 fps

4 Bar code reading for efficient sample control.

Bar code and 2D code (QR code) are easily read. Sample numbers can also be displayed and saved. No more mixed-up samples.



10

NETWORKACCESS

Data sharing through remote access*.

Users (on contract) in remote locations may access and operate the system. Information and data are more widely shared, and medical treatment discussion gets easier. Immediate and accurate observation is now possible.

*For this function, separate contract must be concluded for using the Remote mode. *For communication environment, contact us.



Sample is loaded by the Remote user.




The same image is reviewed both by the user on contract and Remote user. User on contract.



The system is operated both by the user on contract (For focusing, positioning, magnification adjustment, etc.)

5G router



ECLIPSE Ui Display


Sample is removed by the Remote user.



SPECIFICATIONS System


diagram and specifications

DIAGRAM



Keyboard Mouse

DIMENSIONS



SPEC

Model	ECLIPSE Ui US / EU	ECLIPSE Ui US / EU		
Microscopy technique	Diascopic bright-field observa	tion (Koehler illumination)		
Observable specimens	One prepared slide			
	Thickness: 0.9 to 1.7 mm (comprising the thickness of the slide glass,			
	specimen, mounting medium, and cover glass)			
	Microscope slide (ISO8037 compliant)			
	Thickness: 0.9 to 1.2 mm			
	Size: 76 mm x 26 mm			
	· Cover glass (ISO8255-compliant)			
	Thickness: 0.17 mm			
	Size: 18 to 60 mm x 18 to 25 mm			
	* Do not use any other types of microscope slide and cover glass.			
Observable range	Macro-image:	Prepared slide allover (75 mm x 26 mm)		
	Micro-image:	Cover glass allover (60 mm x 25 mm or larger)		
Optical system	Light source:	High-intensity white LED		
		(for macro- and micro-images)		
	Number of fields of view:	22		
	Objectives:	CFI Plan Fluor 4X, CFI Plan Fluor 10X, CFI Plan		
		Fluor 20X, CFI Plan Fluor 40X		
Motorized Function	Nosepiece movement, stage movement (incl. sample loading), objectives focus,			
	aperture stop, sample holder clips open/close, macro/micro observation switching			
Focus drive	Objectives vertical movement system			
	Stroke:	10.3 mm		
	Focusing speed (maximum):	0.7 mm/s or more		
Stage	Stroke:	X; 78 mm, Y; 28 mm		
	Moving velocity (maximum):	X: 85 mm/s, Y: 78 mm/s		
Aperture diaphragm unit	Aperture diameter:	Ø1.2 mm to Ø28.1 mm		
		25% to 100%Supported specimen		
Barcode	2D barcodes:	QR code, Data Matrix code		
	1D barcodes:	CODE-128		
Micro-imaging	X-Y staging:	Half, Full, and Repeat (For Repeat, steps 1 to 6 can		
		be selected for each objective.)		
	Displaymodes:	Live; 1080 x 1080 30 fps		
		High-Quality; 3712 x 3712, 1 fps		
Focus	Contrast AF, manual focus			
Exposure compensation	-2 to +2 EV			
Color adjustment	Contrast, Brightness, Saturation, Hue			
	Point registration*:Registrable up to 12 points			
Marking mode*	Maximum number of spots:	499		
	Spot size:	8px, 16px and 32px selectable		

Measurement mode*	Line segment and length scale display between two points				
Image capture*	Static macro image capture:	1330 x 460			
	Static micro image capture:	Live;	1080 x 1080		
		High-Quality;	3712 x 3712		
	Z-stackimage capture:	Number of images:	1 to 21		
		Photographing interval:	0.5 to 5 µm		
		Setting interval;	0.5 μm		
	Images: Macro (Overall image of observed samp				
	micro (Microscopic appearance)				
Image storage format*	JPEG				
Video capture format*	File format:	MP4			
	Image capturing time:	10 minutes (Maximum)			
	Compression method:	MPEG-4 Video			
	Frame rate:	10 fps			
	Resolution:	1080 x 1080			
Remote operation	Video relay:	AWS, WebRTC			
External interface	LAN:	GbE			
		1000 Mbps (two ports)			
	USB:	USB2.0 cable, Type A			
		480 Mbps (two ports)			
		Mini Display port:	Recommended monitor		
			resolution 1920 x 1200		
Operating system	Windows 10 IoT Enterprise LTS	C 2019			
Main body ratings	Input ratings:	AC100-240 V±10%, 50/6	0 Hz		
	Maximum power consumption: 170 W				
Power cord	• For use in a 100 - 120 VAC region outside Japan:				
	UL-listed detachable cord set, 3-conductor grounding				
	(3-conductor grounding, Type SVT, No. 18 AWG, maximum length 3 m,				
	Plug Type NEMA5-15P, rated at 125 VAC minimum)				
	• For use in a 220-240 VAC region:				
	EU/EN listed detachable cord set, 3-conductor grounding				
	(3-conductor grounding, Type H05VV-F 1 mm2, maximum length 3 m,				
	rated at 250 VAC minimum)				
	• For use in Japan:				
	PSE approved detachable cord set, 3-conductor grounding				
	(3-conductor grounding, Type VCTF 3 x 0.75 mm2, maximum length 3 m,				
	rated at 125 VAC minimum))				

* Not for use in diagnostic procedures.