

# ECLIPSE Ti2-U **Inverted Microscope for ICSI** (I) The second Nikon П Nikon ECLIPSE Ti2 Shedding New Light On MICROSCOPY

## **Superior optical performance and** precision for ICSI

Nikon microscopes have played a significant role from the early days of Intracytoplasmic Sperm Injection (ICSI) and continue to contribute to the advancement of this field through superior optical performance and precision.

The ECLIPSE Ti2-U IVF EU/US provides an ultra-stable platform for performing precision manipulation procedures such as ICSI. The intermediate magnification switching function (1.5X) allows easy switching of observation magnification without changing objectives.



## Nikon Advanced Modulation Contrast (NAMC) observation

The Nikon Advanced Modulation Contrast (NAMC) technique enables high-contrast observation of colorless and transparent samples through plasticware. NAMC produces a pseudo-relief shading effect, ideal for ICSI. The direction of shading can be adjusted through 360 degrees by turning the modulator in an NAMC objective. CFI S Plan Fluor ELWD NAMC objectives are equipped with a modulator ring clamping mechanism to prevent accidental rotation of the modulator ring while rotating the correction collar.



A sperm cell can be clearly observed in the injection pipette Images courtesy of: Derek Keating, B.A. and Gianpiero D. Palermo, MD, PhD., Andrology and Assisted Fertilization: Weill Cornell Medicine Ronald O. Perelman and Claudia Cohen Center for Reproductive Medicine



CFI S Plan Fluor ELWD NAMC 20XC (left), CFI S Plan Fluor ELWD NAMC 40XC (right)

## **Spindle observation**

This polarized light technique enables finely detailed observation of spindle bodies for oocyte assessment and the elimination of spindle damage during sperm injection. Switching between spindle observation and NAMC is fast and easy. The following objectives are compatible with spindle observation: CFI S Plan Fluor ELWD NAMC 20XC and CFI S Plan Fluor ELWD NAMC 40XC.



Spindle bodies (indicated by arrows) observed in different colors depending on their orientation. The color can be changed from blue to red by rotating the condenser module 90 degrees.



TI2-C-SO Spindle observation system



CFI Achromat NAMC 10XF (left) CFI Achromat LWD NAMC 20XF (middle), CFI Achromat LWD NAMC 40XC (right)





## **Specifications**

	Model name	ECLIPSE Ti2-U IVF EU	ECLIPSE Ti2-U IVF US	
Main Body Tube Eyepiece Tube Eyepiece Transmitted Illumination	Optical System	Infinity-corrected CFI60		
	Field Number	22 with C-mount, 25 with F-mount		
	Intermediate Magnification Switching	Manual switching of 1.0x/1.5x		
	Output Port	4 Manual positions Eyepiece 100%, left 100%, right 100%, eyepiece 20%/left 80%		
	Focusing Unit	Manual drive, Coarse/Fine focusing knob, Focusing stroke: 10 mm		
Tube		TC-T-TS Binocular S tube (Field number 22)		
Eyepiece Tube Base		TI2-T-BS S Eyepiece Tube Base Unit		
Eyepiece		CFI 10X Eyepiece (Field number 22)		
Transmitted Illumination	Pillar for transmitted illumination	• TI2-D-PDI Pillar I for Dia Illumination Condenser vertical stroke: 66 mm, Backward tilting up to 25 degrees, with field diagram and refocusing mechanism		
	LED Lamphouse for Dia Illumination	• TI2-D-LHLED LED Lamp House for Dia Illumination: High power LED		
Condenser	Condenser Turret	TC-C-TC Condenser Turret: 7 manual positions		
	Condenser Lens	• TI-C NAMC Condenser Lens (O.D.=44 mm, NA=0.4)		
	Modules for Condenser Turret	TC-C-MN-10X NAMC 10X Module TC-C-MN-20X NAMC 20X Module TC-C-MN-40X NAMC 40X Module		
	Polarizer for NAMC Observation	TC-C-DICPNI NAMC/IMSI Polarizer		
Stage		TC-S-SRI Stage IVF Stroke X: ±57 mm, Stroke Y: ±36.5 mm TI2-S-HL Long Handle	TC-S-SR Stage Stroke X: ±57 mm, Stroke Y: ±36.5 mm TI2-S-HL Long Handle	
Nosepiece		TI2-N-ND DIC Sextuple Nosepiece: 6 manual positions, Simple waterproof structure TI2-N-N Sextuple Nosepiece: 6 manual positions, Simple waterproof structure		
Spindle Observation System		TI2-C-SO Spindle Observation System		
Adjustment Tool		C-T Centering Telescope		
Objectives (NA/W.D.) W.D.=mm		CFI Achromat 4X (0.10/30.0) CFI Achromat NAMC 10XF (0.25/6.20) CFI Achromat LWD NAMC 20XF (0.40/3.10) CFI Achromat LWD NAMC 40XC (0.55/2.70-1.70)	CFI Plan Fluor 4X (0.13/17.20) CFI S Plan Fluor ELWD NAMC 20XC (0.45/8.20-6.90) CFI S Plan Fluor ELWD NAMC 40XC (0.60/3.60-2.80)	
Power Cord and AC Adapter		Power Cord Type BE, 220-240V 4-AC AC Adapter	• Power Cord Type BU,100-120V • I-AC AC Adapter I	
Operating En	vironment	Temperature: 0°C+40°C, Humidity: 60% RH max. (at +40°C	, no condensation), Indoor use only	

#### July 2024 ©2024 NIKON CORPORATION

	TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.	
Monitor images are simulated. Company names and product na	mes appearing in this brochure are their registered trademarks or trademarks.	•

N.B. Export of the products' in this brochure is controlled under the Japanese Foreign Exchange and Foreign Tra Appropriate export procedure shall be required in case of export from Japan. "Products: Hardware and its technical information (including software)





#### NIKON CORPORATION Head office

1-5-20, Nishioi, Shinagawa-ku, Tokyo 140-8601, Japan https://www.healthcare.nikon.com/en/

Manufacturer Japan 471, Nagaodai-cho, Sakae-ku, Yokohama, Kanagawa 244-8533, Japan

#### Nikon Instruments Inc.

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A. phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only) fax: +1-631-547-0299

https://www.microscope.healthcare.nikon.com/

Nikon Europe B.V. (headoffice) Tripolis 100, Burgerweeshuispad 101, 1076 ER Amsterdam, Netherlands phone: +31-20-7099-000 fax: +31-20-7099-298

https://www.microscope.healthcare.nikon.com/en\_EU/

### Nikon Canada Inc.

CANADA phone: +1-905-625-9910 fax: +1-905-602-9953 Nikon France, Succursale de Nikon Europe B.V. FRANCE phone: +33-1-4516-4516 fax: +33-1-4516-4505 Nikon Deutschland, Zweigniederlassung der Nikon Europe B.V.

GERMANY phone: +49-211-9414-888 fax: +49-211-9414-322 **Nikon Italy, Branch of Nikon Europe B.V.** ITALY phone: +39-055-300-9601 fax: +39-055-300-993 Nikon Europe B.V., Amsterdam, Zweigniederlassung Schweiz (Egg/ZH)

SWITZERLAND phone: +41-43-277-2867 fax: +41-43-277-2861 **Nikon UK, Branch of Nikon Europe B.V.** UNITED KINGDOM phone: +44-208-247-1717 fax: +44-208-541-4584 **Nikon Österreich, Zweigniederlassung der Nikon Europe B.V.** AUSTRIA phone: +43-1-972-6111 fax: +43-1-972-6111-40