Nikon

Imaging Software

Elements

Shedding New Light On MICROSCOPY



NIS-Elements is the total imaging solution for your research


Nikon's universal software platform, NIS-Elements, combines powerful image acquisition, analysis, visualization and data sharing tools. With fully customizable user interfaces and software modules, NIS-Elements can serve as a simple interface for photo-documentation and also power complex, conditional workflows with automated imaging and analysis routines.

In addition, utilizing deep learning, a subfield of AI technology, NIS-Elements can perform high-speed image processing and analysis according to the user's specific needs.

Triggering	Confocal	Tracking	
Deep Learning	Large Image Stitching	Deconvolution	Elements
Super-Resolution	Volume Analysis	TIRF	Photostimulation
3D Visualization	Digital I/O	Photodocumentation	Ca ²⁺ Imaging
Time-Lapse	High Content Analysis	Multiphoton	Optogenetics

One software platform for all imaging systems

Nikon believes that having a single software platform for all imaging modalities is vital. NIS-Elements provides the same interface, control, workflow, and terminology whether it's used for widefield, confocal, or super resolution imaging. With one platform to learn, users can easily switch between microscope systems when their applications require different imaging modalities. Imaging results from different Nikon systems can also be easily combined and analyzed to expand your research direction.






Total acquisition-to-analysis solution for high-throughput screening assay

The NIS-Elements HC (High-Content Analysis) provides streamlined high-speed, automated well-plate acquisition, data review, analysis and management for multiple well-plate experiments.


The HC interface simplifies experiment setups using setup wizards. Users can easily define acquisition parameters including well-plate configuration, plate handling, autofocusing, filter switching and detectors.




...

Completely customize to your research

From individual hardware selection and optimization to fine-tuning acquisitions routines and custom multi-channel binary analysis - you are in complete control of tailoring and creating a system built and inspired by your imagination.



LASER LINES	405NM	488NM	
	5%		
		30%	
		GREEN	





Simple customization of advanced analysis automation

The NIS-Elements General Analysis (GA) option enables easy customization of complex analysis or statistical flows such as 3D volume and 4D tracking by simply dragging and dropping analysis templates, ensuring accurate and reliable analyses.




Toward increasing intelligence

Nikon

NIS-Elements allows high-precision, intelligent image processing and analysis using cutting-edge Al technology that employs deep learning. Advanced processes, such as distinguishing samples from noise and extracting cells from unstained samples, which in the past have required a great deal of time and experience, are now carried out automatically, dramatically increasing experiment analysis efficiency and reducing cell damage.

-13823.578

-19927.392



Share your data

NIS Elements

NIS-Elements is designed to get your data "out". There are many options for file and data export to move files, metadata, and analysis results to other formats, other software platforms and even data sharing between programs to leverage other components of your research routines.

Evolves with your research

The software is on the move, always transforming with the demands of research. With NIS-Elements, you can continue to grow your system over time (e.g. upgrade the detector, add additional detectors, change light sources, add a confocal, add high-throughput functionality, etc.).



Nikon's flagship NIS-Elements package

Optimized for advanced research applications, Nikon's flagship software package features fully automated image acquisition, advanced device control and powerful analysis and visualization tools.



Standard research application package

Developed for standard research applications such as analysis and photodocumentation of fluorescent imaging, NIS-Elements BR features up to four-dimensional acquisition and advanced device control capabilities.



Photodocumentation package

Software package for photodocumentation. Include basic measuring and reporting tools.







Confocal imaging package

Dedicated interface for Nikon's confocal and multiphoton systems, providing easy instrument setup and streamlined operation. Incorporates many of the features of NIS-Elements AR for advanced acquisition, image processing, analysis, visualization and data sharing capability.

Confocal resolution enhancement

Higher resolution confocal images can be easily generated with a single click. The software assesses the captured image and automatically determines processing parameters to achieve enhanced resolution.

High Content Analysis

Total acquisition-to-analysis solution for high-content imaging applications. Seamless workflow from microscope and peripheral device control to data analysis and management.

Package Comparison

	CAPTURE	AR	BR	D	С	ER	JOBS/HC
Confocal support					1	1	option
	Time Lapse	1	1	1	1	1	1
	Z-Stack*	\checkmark	✓	1	1	1	1
Nultidio oncional Imagina	Multi Point*	\checkmark	✓	1	1	1	1
Multidimensional Imaging	Multichannel*	1	✓		1	1	1
	4D with Experimental Preview		option				
	6D with Experimental Preview	option			1	1	✓
	AVI Acquisition	\checkmark	\checkmark	1	1	1	1
	JOBS Acquisition	option			option	option	option
	Simultaneous Dual / Triple / Quad Camera	option			option	option	option
Acquisition	Triggered Device Control	option			option	option	option
	DAQ (TTL/ Analog) Control	option	option		option	option	option
	Incubation	option	option		option	option	option
	Volume Contrast	option			option	option	option

* Note: Drivers for third party device control/automation are required.

DISPLAY & PROCESSING		AR	BR	D	С	ER	JOBS/HC
AI	Enhance.ai	Option			Option	Option	Option
	Convert.ai	Option			Option	Option	Option
	Segment.ai	Option			Option	Option	Option
	Denoise.ai	1			1	1	1
	Clarify.ai	Option			Option	1	Option
Image	Annotation	1	1	1	1	1	1
	Image Filters, Morphology	1			1	1	1
	Image Arithmetic	1	1		1	1	1

Image

D/3D

nterac mage mage

1..l+i f

lacro


User m

DISPLAY & PROCESSING		AR	BR	D	C	ER	JOBS/HC
lavers	Multi-dimensional image viewer	1	\checkmark	✓	1	1	1
layers	Binary Layers	1	\checkmark	option	1	✓	✓
DISP image creation ctive display / manipulation	Snapshot	1	\checkmark	1	1	✓	1
	Movie	1	\checkmark	✓	1	✓	1
	Interactive Movie / Volume Rendering	1	\checkmark		1	✓	1
	Tiling (Montage) View	1	\checkmark	✓	1	✓	1
	Max / Min Projections	1	\checkmark		1	✓	1
	Ratio Viewing and Graphing	1			1	✓	1
image creation ctive display /	Plate View, Heat Maps, Sample Labeling	option			option	option	1
	Volume View: 3D ND Crop	1	\checkmark		1	✓	1
	Manual Channel Alignment	1	\checkmark	✓	1	1	1

CAPTURE, DISPLAY & MULTIFUNCTION		AR	BR	D	С	ER	JOBS/HC
	Live Compare	✓	option	option	\checkmark	✓	1
	HDR (High Dynamic Range)	1	option	option	1	✓	1
	EDF / Real Time EDF	option	option	option	option	option	option
functional imaging	2D Large Image Stitching (Free shape)	✓	✓	1	1	✓	1
	3D Large Image Stitching (Free shape)	1	\checkmark		1	✓	1
	FRET/Custom Equation Editor	option			option	option	option
	Deconvolution (2D Real Time/2D/3D)	option			option	1	option
	Macro Creation	1	\checkmark	1	1	✓	1
	Macro Debugger & Variable View	1	option	option	1	1	1
nanagement	Multi-User Environment	1	\checkmark	1	1	1	1

CAPTURE, I	DISPLAY & MULTIFUNCTION	AR	BR	D	С	ER	JOBS/HC
Database	High Content Database	option			option	option	1
Database	Image Database (non HC)	option	option	option	option	option	option
Report	Report Generator	1	1	1	1	1	1

	MEASUREMENT	AR	BR	D	С	ER	JOBS/HC
General measurement Multi-dimentional measurement Tracking Classifier High content	Segmentation	1	1		1	1	1
	Auto-Measurement	1	1	option	1	1	1
	ROI(Region of Interest) Tools & Statistics	1	1		1	1	1
	Time-Measurement	1	option		1	1	1
	Volume Measurement	1			1	1	1
	3D Volume Measurements	option			option	option	option
incusurement	Z profile & 3D EDF Measurements	option	option	option	option	option	option
	Kymograph	1			1	1	1
Tracking	2D/3D Object Tracking	option			option	option	option
	Pixel Classifier	1	1	1	1	1	1
Classifier	Object Classifier-Advanced Segmentation	option			option	option	option
	Colocalization	1			1	1	1
	Live/Dead	option			option	option	1
High content	High Content Intensity Image	option			option	option	1
	Cell Counting	option			option	option	1
	General Analysis (Automated image analysis)	option			option	option	1
Industrial	Grain Sizing, Cast Iron & Filter Analysis	option	option	option	option	option	option



NIS-Elements is not for clinical diagnostic use.

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. May 2023 ©2021-23 NIKON CORPORATION

WARNING TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.

N.B. Export of the products* in this catalog is controlled under the Japanese Foreign Exchange and Foreign Trade Law. Appropriate export procedure shall be required in case of export from Japan. *Products: Hardware and its technical information (including software)

Monitor images are simulated. Company names and product names appearing in this brochure are their registered trademarks or trademarks.

Nikon

NIKON CORPORATION Shinagawa Intercity Tower C, 2-15-3.

Shinagawa Intercity Tower C, 2-15-3, Konan, Minato-ku, Tokyo 108-6290, Japan phone: +81-3-6433-3705 fax: +81-3-6433-3785 https://www.healthcare.nikon.com/en/

Nikon Instruments Inc.

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A. phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only) fax: +1-631-547-0299 https://www.microscope.healthcare.nikon.com/

Nikon Europe B.V.

Stroombaan 14, 1181 VX Amstelveen, The Netherlands phone: +31-20-7099-000 https://www.microscope.healthcare.nikon.com/en_EU/

Nikon Precision (Shanghai) Co., Ltd. CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060 (Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026 (Guangzhou branch) phone: +86-20-3882-0550 fax: +86-20-3882-0580 https://www.nikon-precision.com.cn/

Nikon Canada Inc.

CANADA phone: +1-905-625-9910 fax: +1-905-602-9953 **Nikon France, Succursale de Nikon Europe B.V.** FRANCE phone: +33-1-4516-4516

Nikon Deutschland, Zweigniederlassung der Nikon Europe B.V. GERMANY phone: +49-211-9414-888

Nikon Italy, Branch of Nikon Europe B.V. ITALY phone: +39-055-300-9601

Nikon Europe B.V., Amstelveen, Zweigniederlassung Schweiz (Egg/ZH) SWITZERLAND phone: +41-43-277-2867

Nikon UK, Branch of Nikon Europe B.V. UNITED KINGDOM phone: +44-208-247-1717 ISO 14001 Certified for NIKON CORPORATION



Nikon Österreich, Zweigniederlassung der Nikon Europe B.V. AUSTRIA phone: +43-1-972-6111

Nikon Singapore Pte. Ltd. SINGAPORE phone: +65-6559-3651 fax: +65-6559-3668

Nikon Australia Pty Ltd AUSTRALIA phone: +61-2-8767-6900 Nikon Instruments Korea Co., Ltd.

KOREA phone: +82-2-6288-1900 fax: +82-2-555-4415



En