



The Next Revolution in Microscopy

A Giant Step Forward in Stereo Microscopy

Nikon offers a broad range of stereo microscopes and accessories, including a research stereo microscope system with the world's highest zoom ratio, superb resolution and bright fluorescence imaging. Also features other versatile parallel-optics type models suitable for various applications and Greenough-type models that are user-friendly and affordable.

	SMZ25 SMZ18		SMZ1270/ 1270i	SMZ800N
Optical system	Parallel-optics type			
Zoom ratio	25:1	18:1	12.7:1	8:1
Zooming range	0.63-15.75×	0.75-13.5×	0.63-8×	1-8×
Total magnification*1 (with standard set*2)	3.15-945× (6.3-157.5X)	3.75-810× (7.5-135X)	3.15-480× (6.3-80X)	5-480× (10-80X)
Working distance*3	60mm	60mm	70mm	78mm
Image capture	0	0	0	0
System expandability	0	0	0	0
Embedded use	_	_	0	0

Index

Stereo Microscopes	
• SMZ25, SMZ18 · · · · · · · · · · · · · · · · · · ·	. 4
• SMZ1270/1270i, SMZ800N·····	8
• SMZ745/745T · · · · · · · · · · · · · · · · · · ·	12
• SMZ445/460, SMZ-2·····	13
(CH) M.W.	
Accessories (for SMZ25, SMZ18)	
Base Unit, Focus Unit, Stand/Focus Mount	14
Objectives, Tubes, Nosepiece/Focus Mount Adapter,	
Stage ·····	15
Controller, Darkfield Observation Accessory, Polarizing	
Observation Accessory, Epi-fluorescence Set	16
Fiber Illuminator Set, Coaxial Illuminator,	
Ring LED Illuminator · · · · · · · · · · · · · · · · · · ·	17

Accessories (for SMZ1270/1270i, SMZ800N, SMZ745/745T, SMZ445/460, SMZ-2)

Nosepieces, Tubes, Eye-level Riser,	
Intermediate Tubes · · · · · · · · · · · · · · · · · · ·	19
• Stages, Observation Attachments · · · · · · · · · · · · · · · · · · ·	20
• Illumination Systems · · · · · · · · · · · · · · · · · · ·	21
• Stands	22
Universal Table Stands/Focusing Mounts	23

• Objectives, Auxiliary Objectives · · · · 18

Specifications/System Diagrams 24 • System Diagrams (SMZ25/18) 25

• Specifications (SMZ25/18) · · · · · · · · · · · · · · · · · · ·	25
System Diagrams	
(SMZ1270/1270i/800N, SMZ745/745T) · · · · · · · · · · · · · · · · · · ·	26
Specifications · · · · · · · · · · · · · · · · · · ·	28

SMZ745/SMZ745T	SMZ445/ SMZ460	SMZ-2	
	Greenough type		Optical system
7.5:1	4.4:1 / 4.3:1	5:1	Zoom ratio
0.67-5×	0.8-3.5× / 0.7-3×	0.8-4×	Zooming range
3.35-300× (6.7-50X)	4-70× (8-35X)/ 3.5-60× (7-30X)	4.8-120× (8-40X)	Total magnification*1 (with standard set*2)
115mm	100mm	77.5mm	Working distance*3
○ (SMZ745T)	_	-	Image capture
_	_	-	System expandability
0	0	0	Embedded use

^{*1} Depends on the combination of eyepiece and objective lens *2 With a 10x eyepiece and a 1x objective

^{*3} With a 1x magnification without auxiliary objective

Research Stereo Microscope

SMZ25/SMZ18

Evolutionary stereo microscope

Nikon has developed a stereo microscope that features a large zoom ratio of 25:1, high resolution and exceptional fluorescence transmission capability. These models meet the increasing needs for imaging systems that span spatial scales from single cells to whole organisms.

World's widest zoom range and highest resolution for a stereo microscope

- First stereo microscope to offer a 25:1 zoom range (SMZ25)
- Both eye paths boast numerical apertures (NA) of up to 0.156, using the SHR Plan Apo 1x objective and SMZ25

Automation and digital imaging

- Motorized focus and zoom operation (SMZ25)
- Imaging Software NIS-Elements enables the use of multiple imaging, processing and analysis modalities, including z-stack capture, time-lapse imaging and EDF image generation

Bright, high-contrast fluorescent images

- Fly-eye lens ensures uniform brightness over the entire field of view even at the lowest magnifications
- Breakthroughs in optical design mean significantly improved signal to noise ratio and crystal clear fluorescent images

Easy to use

- User-friendly remote control (SMZ25)
- Easy-to-operate slim LED DIA base with OCC illumination
- · Wide range of illuminators and accessories that accommodate a variety of observation methods





Motorized zoom model with the highest zoom ratio and resolution in Manual zoom model providing advanced optical performance and incredibly bright fluorescence at an attractive price


		r	
Model	SMZ25	SMZ18	
Туре	Motorized zoom	Manual zoom	
Observation	Brightfield/Darkfield/Fluorescence/Simple polarizing		
Zoom ratio	25:1	18:1	
Magnification range	0.63x - 15.75x	0.75x - 13.5x (with 0.75/1/2/3/4/5/6/8/10/12/13.5x click stops)	
Maximum magnification	315x* ¹	270×*1	
Maximum FOV	ø70 mm*²	ø59 mm*²	
Maximum NA of	0.312* ³	0.3*3	

Remarkable resolution and the world's widest zoom range

Dynamic zoom ratio of 25:1 SMZ25







An innovative optical system known as "Perfect Zoom Optics" offers the world's first zoom ratio of 25:1 (zoom range: 0.63x -15.75x). The SMZ25 can seamlessly capture the entire dish while simultaneously delivering microscopic details.

Auto Link Zoom (ALZ) supports seamless viewing at different scales





ALZ automatically adjusts the zoom factor to maintain the same field of view when switching objective lenses. This function enables seamless switching between whole organism imaging at low magnifications and detailed imaging at high magnifications.

Superior resolution never before seen on a stereo microscope SMZ25 SMZ18

Newly developed SHR (Super High Resolution) Plan Apo series objectives offers a resolution of 1100LP/mm (observed value, using SHR Plan Apo 2x at maximum zoom). The 0.5x, 1x, or 1.6x lower magnification objectives deliver a bright field of view and brilliant images with true-to-life colors. All the SHR lenses are parfocal with each other. When used together with the accessory nosepiece,









Comparison of resolution and color aberration by resolution chart

Parallel-optics type


Bright, high-contrast fluorescent images SMZ25 SMZ18

Enhanced brightness and uniform illumination in a low magnification range

The SMZ25 series is the first stereo microscope in the world to use a fly-eye lens on an epi-fluorescence attachment. This ensures bright, uniform illumination even at low magnifications across a large field of view.

Improved S/N ratio and crystal clear fluorescent images thanks to an improved optical system

Nikon has succeeded in improving the signal and reducing noise in fluorescent images by using a short-wavelength, high-transmission Fluor lens. This enables observations of cell division and samples with weak fluorescence, both of which are difficult to observe and record images using conventional stereo microscopes.



Sample images

Fertilized mouse egg Image courtesy of Kazuo Yamagata, Ph.D., Center for Genetic Analysis of Biological Responses, Research Institute for Microbial Diseases. Osaka University



Tg (isl1-GFP) (using SHR Plan Apo 1x at zoom magnification of 6x with SMZ25) Image courtesy of Hisaya Kakinuma, Ph.D., Laboratory for Developmental Gene Regulation, Developmental Brain Science



Automation and digital imaging SMZ25 SMZ18

A wide range of digital imaging capabilities with the Digital Sight series and NIS-Elements imaging software

Easily obtain the information required, such as Z drive position, zoom factor, objective lens, filter cube and LED DIA brightness, by using the Digital Sight series and NIS-Elements together with the microscope.



Detected observation condition/available control ©: Detection and control of observation condition possible O: Detection of observation condition possible				
	SMZ25 + NIS-Elements	SMZ18 + NIS-Elements		
	Motorized focus unit Motorized epi-fluorescence set (control box A)	Manual focus unit Manual epi-fluorescence set (relay box and control box B)		
Zoom magnification	0	0		
Focusing	0	_		
Objective (with nosepiece)	0	0		
Diascopic LED illumination stand (ON/OFF, light intensity control)	0	0		
Fluorescence illuminator (light intensity control)	0	0		
Filter cube	0	0		

For other combinations, please confirm with Nikon.

Note: With NIS-Elements L and F, functions above are not available. Use NIS-Elements D/Br/Ar.

Improved observation efficiency

Easy-to-use OCC illumination SMZ25 SMZ18





The LED DIA Base with built-in OCC illuminator generates minimal heat, consumes little power and has a long life. The illuminator also enhances the contrast of uneven surfaces, such as those of film.



The OCC illuminator can be controlled using a slide lever. Thanks to scales on the slide lever, the user can save and reproduce desired illumination levels. In addition, an OCC plate can be inserted into the illumination unit from the front and rear sides, so images with different shadow direction can be observed.

Comparison images (film) Conventional diascopic illumination OCC illumination

What is OCC illumination?

OCC stands for oblique coherent contrast, a form of oblique lighting method developed by Nikon. Compared to conventional diascopic illumination that illuminates directly from below, OCC illumination applies coherent light to samples in a diagonal direction, adding contrast to colorless and transparent sample structures.

User-friendly remote control SMZ25

The remote control provides easy access to zoom and focus controls and is designed for both right- and left-hand use. The remote control contains an LCD monitor with an adjustable backlight that provides at-a-glance information about zoom factor, objective lens, filter cube and LED DIA brightness.






On-axis imaging for digital images SMZ25

Nosepiece by simply moving the objective lens.



Easily switch between stereo position (stereoscopic view) and mono position (on-axis view) when using the P2-RNI2 Intelligent



Parallel-optics type

Stereo Microscope

SMZ1270/1270i SMZ800N





Incredible sharpness throughout a wide magnification range

These versatile stereo microscopes provide both excellent optical performance, such as high-magnification, high-zoom ratio and high-resolution images, and advanced operability. The expandability of parallel optics makes these models suitable for a wide range of applications.

Highest-in-class zoom ratio

- Highest-in-class zoom ratio of 12.7:1 (0.63 8x) with SMZ1270/1270i
- New WF series objectives optimized for wide viewfield observation at low magnification

High-quality images

 High-level chromatic aberration correction provides sharp images throughout the viewfield.



SMZ1270

Versatile stereo microscope with the highest-in-class zoom ratio

Easy to get results

- Automatically detects magnification data in combination with the digital camera control unit (SMZ1270i only)
- Nosepiece offers both widened magnification range and onaxis imaging
- Eyepiece tubes with various inclination angles and slim-type stands minimize user fatigue during observation

Expandable with a wide range of accessories

 A wide range of accessories are available, including eyepiece tubes and stands that are equal to superior specification stereo microscope models



SMZ1270i

The same as the SMZ1270 but equipped with intelligent functions found in superior models (SMZ1270i with a trinocular tilting tube and nosepiece)



SMZ800N

Affordable model with improved operability and basic performance

Highest-in-class zoom ratio

Wide zoom range

The SMZ1270/1270i offers the highest-in-class zoom ratio of 12.7x (0.63 – 8x). It offers both low-magnification wide viewfield observation of the whole of a 35 mm petri dish* during screening and high-magnification observation of minute cell structures

* with 1x objective at the lowest magnification.



SMZ1270/1270i enables observation of the whole of a 35 mm petri dish.


Wide viewfield of SMZ1270/1270i

Ampulex

0.63x zoom

8x zoom

The SMZ800N comes with a 1-8x zoom range, with higher magnification than conventional models and enables high-resolution observation of 640LP/mm (using ED Plan Apo 2x/WF at maximum zoom).



Wide field objectives

The WF series objectives offer uniformly bright images even at low magnification and wide viewfield observation with SMZ1270/1270i. In addition, a 0.75x objective is now available, expanding the lineup of low magnification objectives.



High-quality images

Apochromat optics have been adopted for the lenses in the SMZ1270/1270i zoom body and semi-apochromat optics in the SMZ800N to achieve high-level chromatic aberration correction. They provide sharp images without blur or color fringe.



Apochromat optics (captured with SMZ1270+Plan Apo 1x/WF)

Conventional optics


Easy to get results

Intelligent function for status readout SMZ1270i

In combination with the imaging software NIS-Elements, the SMZ1270i can detect zoom magnification data. In addition, with the Intelligent Nosepiece P-RNI2 attached, data related to the objective in use is also detected. Calibration data is automatically altered, following changes of magnification, to display the appropriate scale and measurement results on the images.



- Detection of zoom magnification and objective information
- Calibration auto change




 $\label{thm:local_problem} \mbox{Note: With NIS-Elements L and F, functions above are not available. Use NIS-Elements D/Br/Ar.}$

On-axis observation with the nosepiece

The double nosepiece offers easy onaxis imaging, enabling observation of the bottom of holes, accurate simple measurement and extended depth of focus (EDF) imaging without distortion.






Ergonomic design

Eyepiece tubes with a range of inclination angles are available for comfortable observation.

They offer the optimum eyelevel to suit each user. In addition, slim-type plain stands and the LED Diascopic Illumination

Stand easily facilitate the presentation and removal of specimens.







With the LED Diascopic Illumination Stand and Fiber Diascopic Illumination Stand, focus control during observation is possible using the dial in front of the base.

Expandable with a wide range of accessories

In addition to conventional accessories, the level of accessories used with superior models is also available for the SMZ1270/1270i and SMZ800N. These include trinocular tubes and slim-type LED diascopic illumination stands. These allow various microscope configurations to suit numerous routine inspections and a range of research and development applications.








Semiconducto

Brain of adult drosophila excited with GFP Image courtesy of Hokto Kazama, Ph.D., Circuit Mechanisms of Sensory Perception, Brain Science Institute, RIKEN

OCC illumination Diascopic brightfield illumin








With the LED Diascopic Illumination Stand and Fiber Diascopic Illumination Stand, image contrast under OCC illumination can be easily adjusted.

OCC illumination boosts the contrast of transparent sample structures. Hemicentrotus pulcherrimus in two-cell stage

	SMZ1270	SMZ1270i	SMZ800N	
Optical system	Parallel-optics type (zooming type)			
Zoom ratio	12.7 : 1		8:1	
Zoom range	0.63 – 8x (0.63/1/2/3/4/6/8x stops)		1 – 8x (1/2/3/4/6/8x stops)	
Total magnification	3.15 – 480x (depending on eyepiece and ol (with coaxial episcopic illuminator: 15 – 540		5 – 480x (depending on eyepiece and objectives) (coaxial episcopic illuminator: 22.5 – 540x)	
Tubes	Eyepiece inclination: 20° (P-B Binocular Tub	Eyepiece inclination: 20° (P-B Binocular Tube) / 0°-30° (P-TERG100 Trinocular Tilting Tube, P-TERG50 Trinocular Tilting Tube), P-T100 Trinocular Tube		
Eyepieces		C-W10xB (F.N. 22), C-W15x (F.N. 16), C-W20x (F.N. 12.5), C-W30x (F.N. 7)		
Objectives	Plan Apo 0.5x/WF, Plan Apo 0.75x/WF, Plan Apo	o 1x/WF, ED Plan 1.5x/WF, ED Plan 2x/WF	Plan Apo 0.5x/WF, Plan Apo 0.75x/WF, Plan Apo 1x/WF, ED Plan 1.5x/WF, ED Plan 2x/WF, Plan 1x, Achro 0.5x	
Working distance	70 mm (with Plan Apo 1x/WF)		78mm (with Plan 1x)	
Weight (approx.)	9.8 kg (with P-B Binocular Tube + P-DSL32 LED Diascopic Illumination Stand)		6.8 kg (with P-B Binocular Tube + C-PSN Plain Stand)	

Please refer to the system diagram (P. 26-27) for accessory combinations.



Greenough Type Stereo Microscope

SMZ745/745T

Superior 7.5x zoom and 115 mm working distance Trinocular optical head type is also available

- The SMZ745/745T boasts a 7.5x zoom that incorporates the Greenough optical system. The zoom range of 0.67x to 5x provides a broad observation range.
- As well as high zoom ratio and magnification, the SMZ745/745T offers an unrivaled 115 mm working distance.
- The SMZ745T incorporates an optical path switching lever that enables easy switchover between eyepiece and camera. A microscope camera can be attached.



SMZ745T (with Microscope Camera DS-Fi3)

Three "A" design

• Air-tight SMZ745

By making joints air-tight, contamination from dust, oil, water and other contaminants is prevented.

Air-tight construction: JIS Degrees of protection provided by enclosures IPX1

Anti-fungal design SMZ745 SMZ745T

Anti-fungal design developed exclusively by Nikon ensures peace of mind when the microscope is used in environments subject to high heat or humidity.

• Anti-electrostatic SMZ745 SMZ745T

Static electricity built up within the microscope is discharged almost instantly, ensuring higher yields.


Anti-static function: 1000-10V, discharge within 0.2 sec.


1000 1000 0 1 2 (sec)

Specifications		
	SMZ745	SMZ745T
Optical system	Greenough type (zooming type)	Greenough type (zooming type), trinocular tube
Zoom ratio	7.5 : 1	
Zoom range	0.67-5x (with 0.67/1/2/3/4/5x stops)	
Total magnification	3.35-300x (depending on eyepiece and auxiliary objective used)	
Straight tube	-	Built-in C-mount 0.55x magnification lens (F.N. 11), compatible with 2/3 in. or smaller CCD
Tubes	Fixed type Eyepiece inclination: 45° Interpupillary distance adjustment: 52-75 mm	
Eyepieces (with diopter adjustment)	C-W 10xB (F.N. 22), C-W 15x (F.N. 16), C-W 20x (F.N. 12.5), C-W 30x (F.N. 7)	
Auxiliary objectives	G-AL 0.5x (W.D. 211 mm), G-AL 0.7x (W.D. 150 mm), G-AL 1.5x (W.D. 61 mm), G-AL 2x (W.D. 43.5 mm)	
Working distance	115 mm (standard)	
Airtight construction	JIS Degrees of protection provided by enclosures IPX1	-
Weight (approx.)	1.6 kg (body)	1.8 kg (body)

F.N.: Field Number

Dimensions





Greenough Type Stereo Microscope


SMZ445/460


Designed for excellent cost performance

- The SMZ445 has a 45° eyepiece tube inclination, and the SMZ460 has a 60° eyepiece tube inclination, which is ideal for embedded use.
- Compact design with ease-of-use and high optical performance.
- ESD protection guards against electrostatic damage to samples.

	TI I
20	-10
SMZ460 with C-FMBN Focusing Mount BN	
SMZ445 with S-PSN Plain Stand/CN	

Specifications		C-F-SN F-Iail -Stall (U/CN)	
	SMZ445		SMZ460
Optical system		n type (zooming type)	
Zooming ratio	4.4:1		4.3 : 1
Zooming range	0.8-3.5x		0.7–3x
Total magnification	4-70×		3.5-60×
Tube		nation: 45° Interpupillary distance adjustment: 54–75 mm oter adjustable for both eyes	Eyepiece inclination: 60° Interpupillary distance adjustment: 54–75 mm Eyepiece diopter adjustable for both eyes
Eyepieces		DxB (F.N. 21), SM 15xB (F.N. 14), SM 20xB (F.N. 12)	
Auxiliary objectives (option)		0.5x), AL7 (0.7x)	
Working distance	100 mm (stan	nm (standard)	
Weight (approx.)	1.0 kg (body)		1.1 kg (body)





Greenough Type Stereo Microscope


SMZ-2

High-resolution optics ideal for inspection, assembly, and measurement

- Compact design with horizontally positioned zooming ring (rotation: 90°)
- \bullet Eyepiece inclination of 45° for comfortable observation



Specifications	
	SMZ-2
Optical system	Greenough type (zooming type)
Zooming ratio	5:1
Zooming range	0.8–4x
Total magnification	4.8-120x (Depending on eyepiece and auxiliary objective used.)
Tube	Eyepiece inclination: 45° Interpupillary distance adjustment: 56-75 mm
Eyepieces (with diopter adjustment)	SM E10xA (F.N. 23, standard), SM 15xB (F.N. 14), SM 20xB (F.N. 12), C-W30x (F.N. 7)
Auxiliary objectives	AL5 (0.6x), AL7 (0.8x)
Working distance	77.5 mm (with standard configuration)
Weight (approx.)	1.6 kg (body), 1.9 kg (stand)





Wide range of dedicated accessories for SMZ25/SMZ18 for all types of observation

Base Unit, Focus Unit, Stand/Focus Mount

Base Unit

Nikon has improved ease of use by moving the controls to the front of the base, including the brightness adjustment dial and the on/off switch.

Fiber DIA base

The Fiber DIA base features condenser lenses that can be switched between low and high magnifications. Furthermore, the OCC illumination system allows high-contrast illumination.



Slim Bases

The slimmer LED DIA Base and Plain Base help increase efficiency of sample manipulation by bringing the level of the sample closer to



Focus Unit

The focus unit is combined with the base unit. Choose from either a manual or motorized focus unit.



Stand/Focus Mount SMZ18

SMZ18 can be mounted on various compact stands using a focus



SHR Plan Apo Objective Series

The SHR Plan Apo series features higher NA, wider field of view and superior flatness and color aberration correction.

These objective lenses can be seamlessly switched because all magnifications have the same parfocal distance. The bayonet mount design allows lenses to be safely and easily removed.



		SHR Plan Apo 0.5×	SHR Plan Apo 1×	SHR Plan Apo 1.6×	SHR Plan Apo 2×
Maximum	SMZ25	0.078	0.156	0.25	0.312
NA	SMZ18	0.075	0.15	0.24	0.3
Working dis	tance	71 mm	60 mm	30 mm	20 mm
Correction ring		_	_	_	3 mm water depth
Wavelength			380-7	00 nm	

Tubes

3 P2-SHR Plan Apo 1.6x

Choose from two types of tilting trinocular tube and one type of low eyelevel trinocular tube. All tubes have a camera port for seamless integration with the Digital Sight series.

4 P2-SHR Plan Ann 2x





2 P2-TERG50 Trinocular Tilting Tube (eyepiece: port 100:0/50:50)



3 P2-TL100 Trinocular Tube L (eyepiece: port 100:0 / 0:100)

1 P2-TERG100 Trinocular Tilting Tube (eyepiece: port 100:0 / 0:100)

Nosepiece/Focus Mount Adapter

Both single and double nosepieces are available.

1 P2-RNI2 Intelligent Nosepiece



Stage

The stage features an XY stroke of 6x4* inches (150 mm x 100 mm) and can be attached to any of the bases, making it effective for capturing large images when used in combination with imaging software NIS-Elements. A sliding stage and tilting stage are also available. *Limited Y travel with 32 mm column bases



P-SXY64 XY Stage

Remote Control

Nikon offers a remote control unit that can be used to operate the microscope and capture images by hand. A footswitch is also available, allowing the user to operate the microscope and capture images by foot, freeing the hands for sample manipulation.



P2-RC Remote Controller



Darkfield Observation Accessory

Darkfield viewing is possible simply by attaching the darkfield unit to the base.

- 1 P-DF2 LED Darkfield Unit
- Shading cover



Polarizing Observation Accessory

The analyzer is attached to the objective and the polarizer to the base or stand to enable polarized viewing.

1 P2-POL Simple Polarizing Attachment



Epi-fluorescence Set

Motorized Epi-fluorescence Set

The fluorescent turret can be operated using the remote control or imaging software NIS-Elements.

- 1 P2-EFLM Motorized Epi Fluorescence Attachment
- 2 Light shading Plate (comes with Fluorescence Attachment)
- 3 P2-EFL Filter Cube (GFP-B/GFP-L/RFP)
- 4 P2-EFLBF Filter Cube (Brightfield)
- 6 P2-CTLA Control Box
- 6 P2-RC Remote Controller 7 P2-CIA QL1x/0.5x 1/4λ Plate





Combinations with SM725

Manual Epi-fluorescence Set

An easy-to-use manual model for Nikon's highperformance epi-fluorescence attachment.

- 1 P2-EFLI Epi Fluorescence Attachment
- 2 Light shading Plate (comes with Fluorescence Attachment)
- 3 P2-EFL Filter Cube (GFP-B/GFP-L/RFP)
- 4 P2-EFLBF Filter Cube (Brightfield)
- 6 P2-CTLB Control Box
- **6** P2-CIA QL1x/0.5x 1/4 λ Plate





Combinations with SMZ18

Fiber Illuminator Set

Flexible Double Arm Fiber Illumination Set

The direction and angle of illumination can be changed to suit the sample by making adjustments with these double arms. The fiber holder position can also be changed to obtain the optimal position for illuminating samples.

- C-FDF Flexible Double Arm Fiber Illumination Unit
- **2**C-FIDH Fiber Holder
- 3 C-FLED2 LED Light Source for Fiber Illuminator



Combination with SMZ18

Ring Fiber Illumination Set

This ring fiber illumination set features an episcopic illumination unit that effectively captures images (can be used with 1x and 0.5x objective lenses).

- P2-FIR Ring Fiber Illumination Unit
- 2 C-FLED2 LED Light Source for Fiber Illuminator



Combination with SMZ18

Coaxial Illuminator

The coaxial light illuminator makes it possible to view light reflected from the surface of a sample. It is ideal for shooting shadow-less images of thick samples.

- 1 P2-CI Coaxial Epi Illuminator
- 2 C-FLED2 LED Light Source for Fiber Illuminator
- 3 P2-CIA QL1x/0.5x 1/4 λ Plate



Combination with SMZ18

Ring LED Illuminator

Ring LED illuminator is equipped with high-intensity, long-life (20,000 hours) LEDs. The illuminator's dial adjusts the intensity of the white

1 P2-FIRL2 LED Ring Illumination Unit



Combination with SMZ25

A variety of accessories are available for stereoscopic observations



Plan Apo WF series

Objectives

A wide selection with various magnifications and working distances is available, including high-NA, high-resolution and wide-viewfield Plan Apo WF series objectives with superior image flatness and chromatic aberration correction.

SMZ1270/1270i	SMZ800N			
Milean Townson	Wiker with the state of the sta	Nikee and the second se	Niker to to our	Nikos TO FINIS FINIS

4 ED Plan 1.5x/WF

6 ED Plan 2x/WF

- 1 Plan Apo 0.5x/WF
- 2 Plan Apo 0.75x/WF
- 3 Plan Apo 1x/WF

Ob	Objectives dis		Zoom magnification	NA	Actual FOV*1
	0.5x/WF	82	0.63x	0.0095	69.8
	U.3X/ VVF	02	8x	0.0525	5.5
Plan	0.75x/WF	107	0.63x	0.0143	46.6
Apo	U./5X/VVF	107	8x	0.0788	3.7
	1x/WF	70	0.63x	0.0190	34.9
	IX/ VVI	70	8x	0.1050	2.8
	1.5x/WF	44	0.63x	0.0285	23.3
ED	1.58/4//	44	8x	0.1575	1.8
Plan	2x/WF 35	35	0.63x	0.0380	17.5
	ZX/VVF	30	8x	0.2100	1.4

^{*1} With C-W10xB eveniece

SMZ800N







Ob	pjectives	Working distance (mm)	Zoom magnification	NA	Actual FOV*1
Achro	0.5x	189	1x	0.0145	44
ACIIIO	U.5X	109	8x	0.0525	5.5
Dlon	1v	70		0.0290	22
riali	Plan 1x 78		8x	0.1050	2.75

^{*1} With C-W10xB eyepiece


Auxiliary Objectives

Microscopes	Auxiliary objectives	Working distance (mm)
	G-AL 0.5x	211
SMZ745/745T	G-AL 0.7x	150
	G-AL 1.5x	61
	G-AL 2x	43.5

Microscopes Auxiliary objectives		Working distance (mm)
SMZ445/460	AL5 (0.5x)	181
SIVIZ445/460	AL7 (0.7x)	127.5
SMZ-2	AL5 (0.5x)	103
SIVIZ-2	AL7 (0.7x)	95

Nosepieces

Double nosepiece with two-objective switchover. Easy changeover from stereo position (stereoscopic view) to mono position (on-axis view) is possible by simply moving the objective lens to the right.



P-RN2 Nosepiece SMZ1270/1270i SMZ800N

Observations with wider zoom ranges are possible by simply switching between two objectives.



P-RNI2 Intelligent Nosepiece SMZ1270i

Enables easy switchover between two objectives. In combination with the Digital Sight series digital camera, it automatically detects the data of objective in use.



Tubes/Eye-level Riser

SMZ1270/1270i SMZ800N

Various ergonomic tubes with different inclination angles enable suitable eye levels to be selected for observation, even when an intermediate tube or illuminator is attached. Trinocular tubes are also equipped with camera ports.

P-B Binocular Tube

20° inclination angle allows observation without having to lean forward and reduces fatigue during long-time operations.



P-T100 Trinocular Tube

The eyepoint height can be raised 106 mm by tilting the eyepieces 180° up. Optical path switching ratio of eyepiece:camera port is 100:0/0:100



P-TERG100/P-TERG50 Trinocular Tilting Tube

Allows continuous adjustment of the eyepiece inclination from 0° to 30°. Optical path switching ratio of eyepiece:camera port is 100:0/0:100 with P-TERG100 and 100:0/50:50 with P-TERG50.



P-TERG100 Trinocular Tilting Tube

P-IER Eye-level Riser

Increases the eyepoint height 25 mm per riser for a total of 50 mm.



Intermediate Tubes

SMZ1270/1270i SMZ800N

Various intermediate accessories are available that can be inserted between the microscope zooming body and tube.

P-IBSS2 Beam Splitter S2

Using a beam splitter and camera adapter, a digital camera can be attached to the binocular eyepiece tube for imaging. Optical path switching ratio of left eyepiece:right eyepiece:camera port is 100:100:0/100:50:50.

P-THSS Teaching Head

Simultaneous observation of the same viewfield is possible between the eyepiece lenses of both teaching head and microscope, making it ideal for educational purposes. The pointer can indicate target points in the viewfield during observation.

P-IDT Drawing Tube

Drawing sample images is possible by simply tracing observed images that are overlaid on top of drawings within the viewfield. The drawings can be removed from view by using the knob to block the light path.

Stages

Stages allow smooth sample movement in order to change viewfield during observation.

C-SSL Dia-sliding Stage

Used for diascopic observation, this sliding stage can be easily moved in the desired direction simply with a light push. Travel range is within ø38mm.



Can be used with the SMZ25 and SMZ18



Circular Floating Stage 2

Used for episcopic observation. Loaded with a sample, the stage can be easily moved in the desired direction simply with a light push to its edges. Travel range is within ø40mm.

SMZ1270/1270i SMZ800N SMZ745/745T SMZ-2

Can be used with the SMZ25 and SMZ18



C-TRS Tilting Stage

This stage has a nonslip sheet and can be tilted 30° from its horizontal position.

SMZ1270/1270i SMZ800N SMZ745/745T SMZ445/460 SMZ-2

Can be used with the SMZ25 and SMZ18



P-SXY64 XY Stage

The stage features an XY stroke of 150 mm x 65 mm. By attaching AZ100 stage adapters, it can be used for various applications. It can be used with both diascopic and episcopic illuminators.

SMZ745/745T SMZ445/460

SMZ1270/1270i SMZ800N



Observation Attachments

Various observation accessories are available that utilize diascopic and episcopic illuminations. They can be used for samples that are difficult to observe using standard illumination.

P-EFL Epi-fluorescence Attachment

Up to four epi-fluorescence filter cubes can be mounted. The flyeye lens provides bright illumination up to the viewfield periphery.

SMZ1270/1270i SMZ800N



P-DF2 LED Darkfield Unit

Equipped with the while light LED as the light source. Simply placing the unit on the stage enables darkfield observation.

SMZ1270/1270i SMZ745/745T



C-POL Polarizing Attachment

Simple polarizing observation is possible by placing the polarizer on the stage while the analyzer is attached to the tip of the objective lens.

SMZ1270/1270i SMZ800N SMZ745/745T



Illumination Systems

Ring Illuminator

Provides a cone of light from above the sample to the center, minimizing unwanted shadow. Suitable for observation of electronic substrates.



C-FIR Plastic Fiber-optics Ring Illuminator

Illuminator is located away from microscope. It enables bright observation with high-intensity light without damaging sample with its heat.

SMZ1270/1270i SMZ800N SMZ745/745T SMZ445/460 * SMZ-2

* G-OBA60 Adapter is required.

Arm Illuminator/Episcopic Illuminator

The direction and angle of the illumination can be changed with simple adjustments of the









C-FID2 Double Arm Fiber Illuminator

It enables bright observation with high-intensity light without damaging sample with its heat. The direction and angle of illumination can be changed using the flexible arms.

SMZ1270/1270i SMZ800N SMZ-2

SMZ745/745T SMZ445/460

C-FDF Flexible Double Arm **Fiber Illumination Unit**

It enables bright observation with high-intensity light without damaging sample with its heat. The direction and angle of illumination can be

SMZ1270/1270i SMZ800N SMZ745/745T SMZ445/460

changed using the fiber holder.

SMZ1270/1270i SMZ800N SMZ745/745T SMZ445/460

C-LSL2 LED Episcopic Illuminator

In combination with C-PSN Plain Stand/CN and C-PSCN Compact Stand/CN, illumination angle flexibility is possible from the back of the microscope.

By attaching arms, flexible change of direction and angle of illumination is possible.

Coaxial Illuminator

Suitable for brightfield observation for high-reflectance flat surface samples such as polished metals and wafers.



P-CI Coaxial Episcopic Illuminator

Coaxial illuminator for parallel optics-type stereo microscopes. Provides high-intensity illumination for the entire view field.

*1/4 λ plate is required

SMZ1270/1270i SMZ800N

Stands



C-PSN Plain Stand/CN C-PSCN Compact Stand/CN

Offers a comfortable work area and allows easy handling of samples. C-PSCN has a small base that saves desk space.



P-PS32 Plain Stand

Features a slim design with a Ø180 mm stage plate and 160 mm width between the pillar and optical axis to boost working efficiency.



C-LEDS Hybrid LED Stand

Both episcopic and diascopic observations are possible and can be conducted simultaneously. The space-saving built-in illuminator can be switched and adjusted with ease.

Type	<u> Episcopic</u>	Episcopic	Episcopic/Diascopic
Illumination method	_	_	Epi-oblique*, brightfield
Built-in filter	_	_	_
Fine focus knob	_	_	_
Observation V	Vith all objectives, at all zoom ranges	With all objectives, at all zoom ranges	With all objectives, at all zoom ranges

^{*} The illumination area is limited by conditions of use.



C-DS Diascopic Stand S

Features a hand rest for comfortable operation. Used in conjunction with C-DSLU2 LED Unit for Dia Illumination Stand.



P-DSL32 LED Diascopic **Illumination Stand**

The OCC illumination system allows colorless and transparent samples to be observed in high relief. Compact slimtype base enhances operation efficiency.



P-DSF32 Fiber Diascopic **Illumination Stand**

Light source is located away from microscope, enabling bright observation with high-intensity light without damaging sample with its heat.

Туре	Diascopic	Diascopic	Diascopic
Illumination method	Brightfield	Brightfield, OCC**	Brightfield, OCC**
Built-in filter	_	Not required (ø45 mm filter slot provided)	NCB11, ND4/16
Fine focus knob	-	Included	Included
Observation magnification	With all objectives, at all zoom ranges	0.5x objective is compatible with zoom magnifications higher than 1.5x.	0.5x objective is compatible with zoom magnifications higher than 1.5x.
Microscopes F	SMZ1270/1270i SMZ800N SMZ745/745T P-DSL32 and P-DSF32 can be used with the SMZ18.	SMZ445/460	


Universal Table Stands/Focusing Mounts

Universal Table Stands G-US1A/G-US2


These stands are handy in microscopy with large samples not loaded onto the standard stand. The microscope zooming body is mounted to the stand arm via a focusing mount. The G-US1A is a table clamp type (table top thickness: 3 to 53 mm).

- Used in conjunction with the C-FMBN Focusing Mount BN on the SMZ1270/1270i/800N/SMZ745/745T/445/460.
- Used in conjunction with the SM Focusing Mount and the G-USA SM US Adapter on the SMZ-2.
- Cannot be used with the SMZ1270/1270i/800N when intermediate tube is mounted on these models.

G-US1A Universal Table Stand 1



G-US2 Universal Table Stand 2




Universal Table Stand US-3

Not only can it be used for a large sample, but this extremely stable stand also easily accommodates intermediate tubes.

- Used in conjunction with the C-FMAN Focusing Mount AN on the SMZ1270/ 1270i/800N/745/745T/445/460.
- Used in conjunction with the SM Focusing Mount on the SMZ-2.


Universal Table Stand US-3



Specifications				
	Universal Table Stand			
Model	G-US1A	G-US2	US-3	
Vertical cross travel	253mm	254mm	229mm	
Horizontal cross travel	260mm		272mm	
Weight (approx.)	4.4kg	23.0kg	30.5kg	
C-FMAN Focusing Mount AN	-	_	0	
C-FMBN Focusing Mount BN	0		_	
C-FMCN Focusing Mount CN	-		_	
SM Focusing Mount)*	0	


○: Possible * G-USA Adapter is required


G-USA Adapter




Focusing Mounts

Various types of focusing mounts are available depending on use. They are used to incorporate stereo microscope bodies into IC bonders or other devices (SM Focusing Mount is for SMZ-2). These mounts can also be used when attaching microscopes to Universal Table Stands.








Unit: mm

	C-FMAN Focusing Mount AN	C-FMBN Focusing Mount BN	C-FMCN Focusing Mount CN	SM Focusing Mount
Focusing area	40mm	50mm	50mm	40mm
Weight (approx.)	0.6kg	0.8kg	1.6kg	0.6kg
Antistatic function	0	0	_	-
Compatible microscopes	Compatible microscopes SMZ1270/1270i/800N/745/745T/445/460			SMZ-2

** Conditions of use vary depending on objective in use. 22 23


System Diagrams (SMZ25/SMZ18)





Specifications/Dimensions (SMZ25/SMZ18)

	SMZ25	SMZ18		
Zooming body				
Optical system	Parallel-optics type (zooming type), apochromatic optical system			
Zoom	Motorized	Manual		
Zoom ratio	25:1	18:1		
Zoom range	0.63-15.75x	0.75-13.5x (with 0.75/1/2/3/4/5/6/8/10/12/13.5x click stops)		
Aperture diaphragm	Zooming body built-in			
Objectives NA, WD (mm)				
P2-SHR Plan Apo 2x	0.312, 20 (with a correction ring for water 0 to 3 mm in depth)	0.3, 20 (with a correction ring for water 0 to 3 mm in depth)		
P2-SHR Plan Apo 1.6x	0.25, 30	0.24, 30		
P2-SHR Plan Apo 1x	0.156, 60	0.15, 60		
P2-SHR Plan Apo 0.5x	0.078, 71	0.075, 71		
Total Magnification (using C-W10xB eyepieces)	3.15-315x (depending on objective used)	3.75-270x (depending on objective used)		
Eyepieces (F.O.V. mm)	• C-W10xB (22) • C-W 15x (16) • C-W 20x (12.5) • C-W 30x (7)			
Tubes (eyepiece/port)	P2-TERG 100 Trinocular Tilting Tube (100/0 : 0/100) P2-TERG 50 Trinocular Tilting Tube (100/0 : 50/50) Inclination angle : 0-30 degree			
	P2-TL100 Trinocular Tube L (100/0 : 0/100) Inclination angle : 10 degree			
Focusing Unit (stroke from objective's parfocal point)	P2-MFU Motorized Focus Unit (up 96 mm/down 4 mm) P2-FU Focus Unit (up 97 mm/down 5 mm)			
Focus mount adapter/nosepiece	P2-FM Focus Mount Adapter P2-RNI2 Intelligent Nosepiece (2 objectives can be attached)	P2-FM Focus Mount Adapter P2-RNI2 Intelligent Nosepiece (2 objectives can be attached) P2-FMDN Focus Mount (for P-PS32/P-DSL32/P-DSF32 stand		
Bases/stand	P2-PB Plain Base P2-DBL LED Diascopic Illumination Base (OCC P-PS32 Plain Stand (only for SMZ18) P-DSL32 LED Diascopic Illumination Stand (only for SMZ18)			
Stages	P-SXY64 Stage C-SSL DIA Sliding Stage C-TRS Tilting Stage			
Observation methods	Brightfield, Epi Fluorescence, Simple Polarizing (with P2-POL Simple I Oblique Lighting	Polarizing Attachment), Darkfield (with P-DF2 LED Darkfield Unit),		
Epi-fluorescence attachments	4 filter cubes mountable, fly-eye lens built-in			
	P2-EFLM Motorized Epi Fluorescence Attachment P2-EFLI Epi I	Fluorescence Attachment		
	P2-FIRL2 LED Ring Illumination Unit			
Episcopic illuminators	Use with fiber light source • P2-Cl Coaxial Epi Illuminator • P2-FIR Ring Fiber Illumination Unit	C-FDF Flexible Double Arm Fiber Illumination Unit		
Episcopic light source	C-FLED2 LED Light Source for fiber illuminator			
Weight (approx.)	32 kg (Motorized Epi Fluorescence Attachment configuration with Trinocular Tilting Tube, Motorized Focus Unit, Intelligent Nosepiece, LED DIA Base and Objectives 1x and 0.5x)	30 kg (Epi Fluorescence Attachment configuration with Trinocular Tilting Tube, Focus Unit, Intelligent Nosepiece, LED DIA Base a Objectives 1x and 0.5x)		
Power consumption (approx.)	30W (Motorized Epi Fluorescence Attachment configuration with Trinocular Tilting Tube, Motorized Focus Unit, Intelligent Nosepiece and LED DIA Base)	10W (Epi Fluorescence Attachment configuration with Trinocular Tilti Tube, Focus Unit, Intelligent Nosepiece and LED DIA Base)		


Dimensions











System Diagrams (SMZ1270/1270i, SMZ800N, SMZ745/745T)














Specifications

Parallel-optics type					
Model	SMZ25	SMZ18	SMZ1270	SMZ1270i	SMZ800N
Optical system	Parallel-optics type (zooming type)		Parallel-optics type (zooming type)		
Zoom ratio	25:1	18:1	12.7:1		8:1
Zoom range	0.63-15.75x	0.75-13.5x	0.63 – 8x		1 – 8x
Total magnification* (When coaxial episcopic illuminator is attached)	3.15-945x (12.5-472x)	3.75-810x (19-405x)	3.15 – 480x (depending on eyepiece (with coaxial episcopic illuminator: 15		5 – 480x (depending on eyepiece and objectives) (coaxial episcopic illuminator: 22.5 – 540x)
Tubes	P2-TERG 100 Trinocular Tilting Tube, P2-TERG 50 Trinocular Tilting Tube, P2-TL100 Trinocular Tube L		P-B Binocular Tube, P-T100 Trinocular Tube, P-TERG 100 Trinocular Tilting Tube, P-TERG 50 Trinocular Tilting Tube		
Eyepiece inclination	P2-TERG 100/50: 0°-30°, P2-TL100: 10°		P-B: 20° P-T100: 10° P-TERG100/50: 0°-30°		
Interpupillary distance adjustment	P2-TERG 100/50: 50 mm or wider P2-TL100: 50–75mm		P-B: 48–75mm P-TERG100/50: 50 mm or wider		
Eyepieces	C-W10xB (F.N. 22), C-W15x (F.N. 16), C-W20x (F.N. 12.5), C-W30x (F.N. 7) (with diopter adjustment)		C-W10xB (F.N. 22), C-W15x (F.N. 16), C-W20x (F.N. 12.5), C-W30x (F.N. 7) (with diopter adjustment)		
Objectives	P2-SHR Plan Apo 0.5x, P2-SHR Plan Apo 1x, P2-SHR Plan Apo 1.6x, P2-SHR Plan Apo 2x		Plan Apo 0.5x/WF, Plan Apo 0.75x/WF, Plan Apo 1x/WF, ED Plan 1.5x/WF, ED Plan 2x/WF		Plan Apo 0.5x/WF, Plan Apo 0.75x/WF, Plan Apo 1x/WF, ED Plan 1.5x/WF, ED Plan 2x/WF, Plan 1x, Achro 0.5x
Working distance (with standard configuration or 1x objective)	60 mm		70 mm 78 mm		78 mm
Weight (approx.)	32 kg (motorized Epi Fluorescence Attachment configuration)	10 kg (with Plain Stand and Ring LED set)	9.8 kg (with Binocular Tube + LED Diascopic Illumination Stand)	11.9 kg (with Trinocular Tilting Tube + LED Diascopic Illumination Stand)	6.8 kg (with Binocular Tube + Plain Stand)

^{*} Depending on eyepiece and objective used

Greenough type						
Model	SMZ745/745T	SMZ445	SMZ460 SMZ-2			
Optical system	Greenough type (zooming type) Trinocular Tube (SMZ745T)	Greenough type (zooming type)	Greenough type (zooming type)			
Zoom ratio	7.5 : 1	4.4 : 1	4.3 : 1	5:1		
Zoom range	0.67–5x	0.8–3.5x	0.7–3x	0.8–4x		
Total magnification*	3.35–300x	4–70x	3.5–60x	4.8–120x		
Tubes	Fixed (binocular tube: SMZ745, trinocular tube: SMZ745T)	Fixed	Fixed			
Eyepiece inclination	45°	45°	60°	45°		
Interpupillary distance adjustment	52–75mm	54–75mm	54–75mm	56–75mm		
Eyepieces	C-W10xB (F.N. 22), C-W15x (F.N. 16), C-W20x (F.N. 12.5), C-W30x (F.N. 7) (with diopter adjustment)	SM 10xB (F.N. 21), SM 15xB (F.N. 14), SM 20xB (F.N. 12)	SM 10xB (F.N. 21), SM 15xB (F.N. 14), SM 20xB (F.N. 12)	SM E10xA (F.N. 23, standard), SM 15xB (F.N. 14), SM 20xB (F.N. 12), C-W30x (F.N. 7)		
Objectives	-	_	-	0.8-4x		
Auxiliary objectives	G-AL 0.5x (W.D. 211mm), 0.7x (W.D. 150mm), 1.5x (W.D. 61mm), 2x (W.D. 43.5mm)	SM-AL 0.5x, 0.7x	SM-AL 0.5x (W.D. 181mm), 0.7x (W.D. 127.5mm)	AL5 (0.6x, W.D. 103mm), AL7 (0.8x, W.D. 95mm)		
Working distance (with standard configuration or 1x objective)	115mm	100mm	100mm	77.5mm		
Weight (approx.)	1.6kg (SMZ745 body) 1.8kg (SMZ745T body)	1.0kg (body)	1.1kg (body)	1.6kg (body), 1.9kg (Stand)		

^{*} Depending on eyepiece and objective used

Related Products

Camera Head

Microscope Camera

Microscope Camera

Digital Sight 1000 Digital Sight 10









Monochrome Microscope Camera

Digital Sight 50M DS-Fi3





Microscope Camera







Elements

Imaging Software Enables a wide range of advanced digital imaging capabilities using a desktop PC and tablet PC.














Simply installing NIS-Elements L on a tablet PC enables setting and control of microscope cameras, live image display, and image acquisition.

Note: NIS-Elements L is available only for

Digital Sight 1000, DS-Fi3, and Digital Sight 10.


EDF (Extended Depth of Focus)

Captures multiple high-resolution images at different focal depths to create a single extended depth of focus image or quasi-3D image.



Multichannel (multicolor)

Multiple fluorescent channels can be captured in conjunction with other imaging methods, such as OCC or brightfield.



Overlapping image with all colors

Individual cells resolved in a live drosophila embryo expressing GFP

(Using SHR Plan Apo 2x at zoom magnification of 8x with SMZ25) Image courtesy of Max V. Staller, Ph.D., Clarissa Scholes, and Angela DePace, Ph.D.,

A wide variety of tools

NIS-Elements L enables the conducting of simple measurements on images, with input of lines and comments. These can also be written onto and saved with the image, and measurement data can be output.





Graticule/scale function							
CrosshairsSimple crosshairsCircle	☐ Grid☐ Horizontal scale☐ Vertical scale						


Optical Systems


Parallel-optics type (zooming type)

This system has a parallel optical path into which various intermediate tubes, including a beam splitter, coaxial episcopic illuminator, epi-fluorescence attachment, teaching head, drawing tube and eye-level riser, can be inserted.

Greenough type (zooming type)

Allows a compact body that is suited for incorporation into other devices.





N.B. Export of the products* in this catalog is controlled under the Japanese Foreign Exchange and Foreign Trade Law. Appropriate export procedure shall be required in case of export from Japan.

*Products: Hardware and its technical information (including software)

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. December 2023 ©2007-23 NIKON CORPORATION



WARNING

TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.



NIKON CORPORATION

Shinagawa Intercity Tower C, 2-15-3, Konan, Minato-ku, Tokyo 108-6290, Japan Healthcare Business Unit, phone: +81-3-6433-3705 fax: +81-3-6433-3785

https://www.healthcare.nikon.com/

Industrial Metrology Business Unit, phone: +81-3-6433-3701 fax: +81-3-6433-3784 https://industry.nikon.com/

NIKON INSTRUMENTS INC.

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A. phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only) fax: +1-631-547-0306

https://www.microscope.healthcare.nikon.com/

NIKON METROLOGY, INC.

12701 Grand River Avenue, Brighton, MI 48116 U.S.A. phone: +1-810-220-4360 fax: +1-810-220-4300 **E-mail: Sales.US-NM@nikon.com**

https://industry.nikon.com/en-us/ NIKON EUROPE B.V.

Stroombaan 14, 1181 VX Amstelveen, The Netherlands phone: +31-20-7099-000

https://www.microscope.healthcare.nikon.com/en_EU/

NIKON METROLOGY EUROPE NV

Interleuvenlaan 86 B-3001 Leuven, Belgium phone: +32-16-74-01-00 fax: +32-16-74-01-03 E-mail: Sales.Europe.NM@nikon.com

https://industry.nikon.com/en-gb/

Nikon Precision (Shanghai) Co., Ltd. CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060 (Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026 (Guangzhou branch) phone: +86-20-3882-0550 fax: +86-20-3882-0580 https://www.nikon-precision.com.cn/

NIKON INSTRUMENTS KOREA CO., LTD.

KOREA phone: +82-2-2186-8400 fax: +82-2-555-4415 https://www.microscope.healthcare.nikon.com/ko_KR/

NIKON SINGAPORE PTE. LTD.

SINGAPORE phone: +65-6559-3651 fax: +65-6559-3668

NIKON AUSTRALIA PTY LTD NIKON INDIA PVT. LTD.

NIKON METROLOGY - MÉXICO MEXICO phone: +52-442-688-5067 E-mail: Sales.NM-MX@nikon.com

NIKON CANADA INC. CANADA phone: +1-905-625-9910 fax: +1-905-602-9953

NIKON ITALY, BRANCH OF NIKON EUROPE B.V.

NIKON EUROPE B.V., AMSTELVEEN,

ZWEIGNIEDERLASSUNG SCHWEIZ (EGG/ZH) SWITZERLAND phone: +41-43-277-2867

NIKON ÖSTERREICH, ZWEIGNIEDERLASSUNG DER

NIKON EUROPE B.V. AUSTRIA phone: +43-1-972-6111

NIKON UK LTD.

UNITED KINGDOM: phone: +44-208-247-1717

ISO 14001 Certified for NIKON CORPORATION

ISO 9001 Certified for NIKON CORPORATION Industrial Metrology Business Unit

NIKON METROLOGY UK LTD.

UNITED KINGDOM: phone: +44-1332-811-349 fax: +44-1332-639-881 E-mail: Sales.UK.NM@nikon.com

NIKON FRANCE, SUCCURSALE DE NIKON EUROPE B.V. FRANCE phone: +33-1-4516-4516

NIKON METROLOGY SARL

FRANCE phone: +33-1-60-86-09-76 fax: +33-1-60-86-57-35

E-mail: Sales.France.NM@nikon.com

NIKON DEUTSCHLAND, ZWEIGNIEDERLASSUNG DER NIKON EUROPE B.V.

GERMANY phone: +49-211-9414-888

NIKON METROLOGY GMBH

GERMANY phone: +49-211-45-44-69-51
E-mail: Sales.Germany.NM@nikon.com