

# Multi-omics spatial analysis of colon cancer tissue reveals emergence of an immunosuppressive tumor maintenance system



Gaurav Joshi<sup>1</sup>, Nicholas Sciascia<sup>1</sup>, Michael Yang<sup>1</sup>, Brian Falkenstein<sup>2</sup>, Raymond Yan<sup>2</sup>, Shannon Quinn<sup>2</sup>, Brenna Fearey<sup>1</sup>, Junya Yoshioka<sup>1</sup>, Akif Burak Tosun<sup>2</sup>, S. Chakra Chennubhotla<sup>2</sup>, Filippo Pullara<sup>2</sup>, Fumiki Yanagawa<sup>1</sup>  
Nikon Instruments, Inc., Boston, MA<sup>1</sup>, PredxBio Inc, Pittsburgh, PA<sup>2</sup>

## Abstract

Using a combination of modalities such as sequential multiplex immunofluorescence (mIF) and spatial transcriptomics (STx) on colon cancer sample we investigated spatial distribution, phenotype, function, and gene expression profile within the colon cancer and healthy tissue. Immunosuppressive cells, Treg and M2 macrophages were in the TME. Presence of other immune cells was not sufficient to prevent crypt hyperproliferation. While immunosuppressive phenotype was observed using antibodies, a strong presence of immunosuppressive genes was not found in transcriptomic study. Analysis of this multi-omics dataset showed hypoxia, a known factor to drive immunosuppression and COL6A2 gene expression to have strong association in driving tumor progression.

## Multi-omics Spatial Imaging and Analysis



## Antibody Panel

Checkpoint	Structural	Stromal	T cells	Autophagy & Hypoxia
PD-1	E-Cadherin	$\alpha$ -SMA	CD3	CA-9
PD-L1	PanCK	Podoplanin	CD4	LC3B
IDO-1			CD8	
VISTA	NK Cell	Neutrophil	FOXP3	Endothelial
ICOS	CD56	CD66b	Treg	CD31
CTLA-4	B Cell	CD177	TCR delta	
TIM-3	CD20		Granzyme B	

5 um FFPE sections were prepared from stage IIA colon cancer sample of a 56-year-old individual exhibiting moderate to poorly differentiated adenocarcinoma of T3, N0 staging. The tumor was microsatellite stable (MSS) and mismatch repair proficient (pMMR). There was low tumor mutational burden.


SpaceIQ™, PredxBio's explainable AI-driven spatial multi-omics platform that provides a rigorous framework to interrogate the TME, is utilized to analyze data. By leveraging probabilistic modeling and unbiased cell typing, SpaceIQ™ uncovered transition cell states, cell-cell interactions, and mechanistic pathway networks within the tumor – yielding a deeper understanding of tumor heterogeneity and therapy resistance.

## References


- Li X, Li Z, Gu S, Zhao X.. A pan-cancer analysis of collagen VI family on prognosis, tumor microenvironment, and its potential therapeutic effect BMC Bioinformatics. 2022 Sep 27;23(1):390 PMID: 36167487

**NIKON BIOIMAGING LAB**  
YOUR SPATIAL BIOLOGY CRO PARTNER  
Contact: gaurav.joshi@nikon.com  
**DOWNLOAD POSTER**

## Results




**FIGURE 1 - Multiplex IF and Visium HD on adjacent sections from colon cancer sample show hyperproliferated crypts, diversity of TME, and gene expression results for various mediators of tumor progression:** E-Cadherin as a marker of epithelial cells is labelled in Cyan highlighting the crypts. Various immune cells and stromal markers surround the crypt in TME. With an exception of CD177+ neutrophils, none of the immune markers used in this study invade the tumor. Stromal proteins such as  $\alpha$ -SMA and Podoplanin surround the hyperproliferated crypt. For the Visium HD data, co-expression matrix of various genes in relation to E-Cadherin shows their distribution in TME.



**FIGURE 2 – Multiplex IF and Visium HD on adjacent sections from a healthy colon tissue show variation in localization and presence of proteins compared to diseased tissue:**  $\alpha$ -SMA+ cells are present around the crypt, at a lower density and further from crypt, while Podoplanin is absent. M2 macrophages and NK cells were present in lamina propria.

## Multi-omics Spatial Analysis with PredxBio's SpaceIQ™ Platform



**FIGURE 3 - SpaceIQ platform identifies COL6A2 gene and hypoxia as key mediators in disease progression:** A pointwise mutual information (PMI) for probabilistic cell-cell association was derived analyzing various features obtained from H&E cell typing. Overlaying the PMI on mIF or Visium HD data allowed us to identify COL6A2 gene and CA9 protein involved in hypoxia as key mediators

## Conclusions

- Comparison between stage IIA, pMMR colon cancer and healthy colon tissue, show structural changes with immune cells and fibroblasts enriched TME. These give an insight in to tumor progression, it's potential to metastasize and help find ways to mitigate it.
- CAFs interact with immune cells to create an immunosuppressive environment promoting cancer growth. Drugs targeting CAFs to inhibit/reprogram them are an active area of research.
- SpaceIQ platform identified CA9, hypoxia marker and COL6A2 gene for having strong feature associations. COL6A9 has been shown to contribute towards tumor progression by vascular remodeling from ECM, promoting drug resistance<sup>1</sup>.

**PredxBio**  
Reimagining Cancer Care with AI-Driven Spatial Biomarkers  
Contact: fillipo@predxbio.com  
**SpaceIQ™ platform**