



Mapping the tumor landscape: Spatial Multi-omics analysis reveals immune and stromal heterogeneity across early and late-stage NSCLC



Gaurav N Joshi¹, Nicholas Sciascia¹, Michael Yang¹, Junya Yoshioka¹, Fumiki Yanagawa¹
Nikon Instruments Inc., Boston, MA¹

Abstract

Using a combination of modalities such as spatial transcriptomics (STx) and sequential multiplex immunofluorescence (mIF) on different stages of non-small cell lung cancer samples we investigated spatial distribution, phenotype, function, and gene expression profile. In spite of immune cell presence stage IIB tumor has evolved and adapted to thrive under pressure with upregulation of PD-L1 on epithelial cells leading to immune evasion. All three tumor types were found to have M2 macrophage population which is prognostically undesirable. The tertiary lymphoid structures had different cell composition between stages IB and IIB tumors, whereas the stage IIIB tumor was devoid of TLS. Unbiased differential gene expression analysis across three tumor types show unique cluster for stage IIIB tumor.

Multi-omics Spatial Imaging and Analysis



10x Visium CytAssist (Gene Expression) + **Lunaphore COMET (Multiplex IF)** + **Data Analysis HORIZON™**

Antibody Panel

Checkpoint	Structural	Stromal	T cells	Other
PD-1	E-Cadherin	α-SMA	CD3	CD11b
PD-L1	PanCK		CD4	Granzyme B
IDO-1		Cell Prolif	CD8	
VISTA	NK Cell	Ki67	FOXP3	Endothelial
LAG-3	CD56			CD31
	B Cell	Macrophage		Immune
	CD20	CD68 CD206 (M2) CD163 (M2)		CD45
				HLA-DR

Non-small Cell Lung Cancer sample distribution

Stage	Pathology	TNM	Grade
I B	Adenocarcinoma	T2N0MX	G1
II B	Adenocarcinoma	T2aN1M0	G1
III B	Adenocarcinoma with micropapillary features	T3N2	G2





NIKON BIOIMAGING LAB

YOUR SPATIAL BIOLOGY CRO PARTNER

Contact: gaurav.joshi@nikon.com


← **DOWNLOAD POSTER**

Results





Top 10 gene expression for the unique green cluster in Stage IIIB NSCLC

Gene	Prognostic value
ANK1	Associated with tumor cell resilience
SCGB3A2	Maintain airway epithelial lineage and modulates inflammatory signaling
FMO5	Enables tumor cells to tolerate ROS and inflammatory oxidative stress. Reflects chronic immune pressure (NK/T-cell activity)
SDR16C5	Supports metabolic plasticity during stress
PCSK1N	Immune modulation
SYBU	Vesicular trafficking, potentially enabling PD-L1 expression
SPTB	Retains lineage identity often seen in immunogenic IFN-γ responsive tumor niches
SFTA2	Chronic inflammatory exposure
PGC	Tumor adaptation
REEP6	Membrane turnover and protein transport, likely facilitating PD-L1 expression



	Stage IB	Stage IIIB
Immunogenicity and Antigen presentation	High considering PanCK+ HLA-DR+ cells	Low considering PanCK+ HLA-DR-cells
TME phenotype	Likely hot	Immune evasion
CD4+ T cells	Antigen recognition and regulatory feedback to control response due to presence of CD4+FOXP3+ population	
CD8+ T cells	Present	
NK cells	Present but do not express Granzyme B	A subpopulation actively express Granzyme B
PD-L1 expression	Absent from epithelial cells	Present on epithelial cells to avoid cytotoxic response

Tumor progression score (TPS) was calculated for Stage IIIB sample on PanCK+ cells expressing PD-L1 . The score was 70%



Conclusions

- A combination of tools, identify unique cell states with clustering helping to identify closeness between different cells which vary between different NSCLC stages.
- We are currently leveraging the SpatialX platform to perform more in depth spatial gene expression analyses to study how the TME changes through NSCLC progression and any changes in gene expression pattern at different intratumoral niches and it's relationship to protein expression by overlaying Visium and mIF datasets.